



### Thermo Scientific Orion VERSA STAR

Advanced Electrochemistry Meter

Visible indicators to quickly show how each of the four channels are assigned

Calibration information with buffer points and electrode slope from last calibration

Function keys that update depending on your screen and shortcut keys on the keypad make this meter easy for anyone to use



Four channels to mix and match in any combination as you need

The large, color graphic display is easy on the eyes and can be seen from any angle

Display information can be customized to show only the information you need

More BIG surprises just for you! Find out on the following pages.

This versatile meter has a large, color graphic LCD display. The informative (and customizable) display shows readings, connections, time, date, sample ID, user ID and calibration information.

- Select from five modules to customize the four channels your way it's versatile for you to configure how you need it for now or for later
- Anyone can use this meter with its plain language screen prompts, function keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Memory for up to 10 methods per channel allow for customized procedures to differentiate between tests and/or users

- Data transfer and updating meter software is easy with two USB and RS232 ports
- Includes data analysis software for transferring stored results or displaying readings in real-time
- Mix samples without a stir plate with direct control for two Orion Star stirrer probes (096019, sold separately)
- Included electrode arm and newlydesigned probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter
- IP54-rated housing handles splashes
- 3 year meter replacement warranty



## A variety of modules that you can place in any of the four channels.

### pH/Temperature Module Allows:

- Up to 6 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Fixing of calibration errors without a complete recalibration with calibration editing

### • Use the pH/ISE/Temperature Module for:

- All the benefits of the pH module
- Fixing pH/ISE calibration errors without a complete recalibration with calibration editing
- Advanced ISE benefits such as auto-blank, timed end point, linear point-to-point, non-linear selectable auto-blank, low-concentration range stability, ISE ID and incremental techniques

### • Conductivity/Temperature Module Offers:

- Conductivity, TDS, salinity or resistivity capabilities
- Selectable reading reference temperatures of 5, 10, 15, 20 or 25 °C with linear, non-linear, nLFu or USP/EP curve options for conductivity
- Choice of practical salinity or natural sea water curves for salinity readings
- Linear or ISO/EN 27888 curves for TDS

### • Dissolved Oxygen/RDO/Temperature Module:

- Takes the worry out of DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- Compensates automatically for temperature, salinity with manual input and barometric pressure (automatic or manual)
- Allows for calibration using water-saturated air, air-saturated water, a Winkler titration or zero point
- Calculates BOD/DOUR/SOUR readings when combined with the complimentary data analysis software

### Module Featuring LogR Technology:

- Includes all of the features of the pH module with the added benefits of:
- The most advanced electrode diagnostics
- No temperature sensor needed to get temperaturecompensated pH readings







### Thermo Scientific Orion VERSA STAR 90

pH/ISE/Conductivity/RDO/Dissolved Oxygen Benchtop Meters

Looking for the top-of-the-line, best meter you can find? Looking for something for everything and meeting all those testing needs? Whether it's one sample that needs a variety of testing or a variety of samples needing one test, the Orion VERSA STAR 90 and 91 advanced multiparameter meters are ready for the job. These meters come fully-loaded and ready to take on multiple measurements.

### **Highlighted Features Include:**

- Display customization to show only the information you need
- Plain language screen prompts, function keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese make it easy for anyone to use this meter
- Non-volatile memory that holds up to 2000 data points with time and date stamp
- Stores up to 10 methods per channel allow for customized procedures to differentiate between tests and/or users
- Two USB and a RS232 port makes for easy data transfer and updating meter software
- Complimentary data analysis software for easy transferring of stored results, real-time viewing and calculation of BOD, DOUR and SOUR results
- Option to control sample mixing directly from the meter without a separate stir plate when using with the Orion Star stirrer probe
- Increase your productivity two stirrer probe ports allow for simultaneous stirring and control for two separate samples at the same time
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Working with almost every AC power source through the included universal power adapter
- IP54-rated housing handles splashes and cleans up easily
- 3 year meter replacement warranty

### **Product Specifications – Orion VERSA STAR Meter**

|                                   |                       | CE   TUV   EPA   |
|-----------------------------------|-----------------------|--|
| Maximum Number of Modules         |                       | 4  |
| Methods                           |                       | Up to 10   |
| Datalogging                       | Number of<br>Datalogs | 2000 with time and date stamp                            |
|                                   | Log Function          | Manual, AUTO-READ™, timed interval, single timed reading |
| Communication<br>Ports            |                       | One RS232, two USB ports (one client, one host)          |
| Orion Star Stirrer<br>Probe Ports |                       | 2  |
| Power                             |                       | Universal, AC power adapter,<br>90-260 VAC, 50-60 Hz     |

### Product Specifications - pH/ISE Module (VSTAR-ISE)

|             |                           | CE   TUV   EPA  |
|-------------|---------------------------|---|
| рН          | Range                     | -2.000 to 20.000  |
|             | Resolution                | 0.1, 0.01, 0.001  |
|             | Relative<br>Accuracy      | ±0.002  |
|             | Calibration<br>Points     | Up to 6   |
|             | Calibration<br>Editing    | Yes   |
|             | Adjustable<br>ISO Point   | Yes   |
|             | Range                     | ±2000.0 mV  |
| mV/Rel mV/  | Resolution                | 0.1   |
| ORP         | Relative<br>Accuracy      | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|             | E <sub>H</sub> ORP Mode   | Yes   |
|             | Range                     | 0 to 19999  |
|             | Resolution                | Up to 4 significant digits  |
|             | Relative<br>Accuracy      | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|             | Units                     | ppm, M, mg/L, %, ppb, none  |
|             | Calibration Points        | Up to 6   |
| ISE         | Calibration<br>Editing    | Yes   |
| IJL         | Adjustable<br>ISO Point   | Yes   |
|             | Calibration<br>Features   | Timed end point, linear point to point,<br>non-linear selectable auto-blank, low<br>concentration range stability     |
|             | Incremental<br>Techniques | Single known addition & subtraction,<br>double known addition & subtraction,<br>single analate addition & subtraction |
|             | Advanced<br>Features      | Displays electrode type, temperature compensation   |
|             | Range                     | -5 to 105 °C, 23 to 221 °F  |
| Temperature | Resolution                | 0.1   |
|             | Relative<br>Accuracy      | ±0.1  |
|             | Offset<br>Calibration     | 1 point   |
| Inputs      | pH/ISE                    | BNC, reference pin  |
| iiiputs     | ATC                       | 8-pin mini-DIN  |

### Product Specifications - Conductivity Module (VSTAR-CND)

|              |                              | CE   TUV   EPA                                     |
|--------------|------------------------------|--|
| Conductivity | Range                        | 0.001 µS to 3000 mS                                |
|              | Resolution                   | 0.001 µS minimum;<br>4 significant figures minimum |
|              | Relative<br>Accuracy         | 0.5 % reading ±1 digit                             |
|              | Reference<br>Temperature     | 5, 10, 15, 20, 25 °C                               |
|              | Temperature<br>Compensation  | Linear (0 to 10.0 %/°C), nLn, nLFu, USP/EP         |
|              | Compatible<br>Cell Constants | 0.001 to 199.9                                     |
|              | Calibration<br>Points        | Up to 6 points                                     |
|              | Calibration<br>Editing       | Yes  |



| Resistivity | Range                 | 2 ohm to 100 meg-ohm                         |
|-------------|-----------------------|--|
|             | Resolution            | 2 ohms-cm                                    |
|             | Relative<br>Accuracy  | 0.5 % reading ±1 digit                       |
|             | Туре                  | Practical salinity or natural sea water      |
| Salinity    | Range                 | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water |
| Samily      | Resolution            | 0.01   |
|             | Relative<br>Accuracy  | ±0.1   |
|             | Range                 | 0 to 200 ppt                                 |
|             | Resolution            | 4 significant digits                         |
| TDS         | Relative<br>Accuracy  | ±0.5 % reading ±1 digit                      |
|             | TDS Factor<br>Range   | Linear 0.01 to 10.00                         |
|             | Range                 | -5 to 105 °C, 23 to 221 °F                   |
| Temperature | Resolution            | 0.1  |
|             | Relative<br>Accuracy  | ±0.1   |
|             | Offset<br>Calibration | 1 point                                      |
| Input       |                       | 8-pin mini-DIN                               |

### Product Specifications - RDO/DO Module (VSTAR-RD)

| Froduct Sp          | ecilications                                      | - NDO/DO WIOGUIE (VSTAN-ND)   |  |
|---------------------|---|---|--|
|                     |   | CE   TUV   EPA  |  |
|                     | Concentration                                     | - Polarographic   |  |
|                     | Range   | 0 to 90 mg/L  |  |
|                     | Resolution  | 0.01, 0.1   |  |
|                     | Relative<br>Accuracy                              | ±0.2  |  |
|                     | % Saturation -                                    | Polarographic   |  |
|                     | Range   | 0 to 600  |  |
|                     | Resolution  | 0.1, 1  |  |
|                     | Relative<br>Accuracy                              | ±2 %  |  |
|                     | Concentration                                     | - RDO   |  |
|                     | Range   | 0 to 50 mg/L  |  |
|                     | Resolution  | 0.01, 0.1   |  |
| Dissolved<br>Oxygen | Relative<br>Accuracy                              | ±0.1 mg/L up to 8 mg/L;<br>±0.2 mg/L from 8 to 20 mg/L;<br>10 % of reading, 20 to 50 mg/L |  |
|                     | % Saturation - RDO                                |   |  |
|                     | Range   | 0 to 500  |  |
|                     | Resolution  | 0.1, 1  |  |
|                     | Relative<br>Accuracy                              | ±2 %  |  |
|                     | Automatic<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg  |  |
|                     | Salinity<br>Factor<br>Correction                  | Automatic when used with conductivity module and probe or manual; 0 to 45 ppt             |  |
|                     | Calibration<br>Features                           | Water saturated air, air saturated water, manual (Winkler) and zero point                 |  |
|                     | Probe Type  | Polarographic or RDO  |  |
|                     | Range   | 0 to 50 °C, 32 to 122 °F  |  |
|                     | Resolution  | 0.1   |  |
| Temperature         | Relative<br>Accuracy                              | ±0.1  |  |
|                     | Offset<br>Calibration                             | 1 point   |  |
| Input               |   | 9-pin mini-DIN  |  |
|                     |   |   |  |

### **Meter Ordering Information**

Description

Cat. No.

Orion VERSA STAR pH/ISE/Conductivity/RDO/Dissolved Oxygen Benchtop Multiparameter Meter Set

Thermo Scientific Orion Laboratory Products Catalog | Sourcebook

рН 📉

ISE 📉

ORP \_

Cond

RDO

DO **T** 





### Thermo Scientific Orion VERSA STAR 40

pH/ISE Benchtop Meters

Measuring pH/ISE has never looked so good or been so easy. The large color display is easy on to read and easy on the eyes. Prompts make accuracy and repeatability easy for even the most complicated testing. If you need to add more channels for more or obtain more parameters from the same, simply remove the protective individual channel cover and add the new module.

### **Features and Benefits**

- Versatility is the name of the game
  - Change out modules to create the meter you need at any time
  - Customize the display to show only the information you need
- Up to 6 point pH calibration with automatic recognition for USA/NIST and DIN buffers and electrode identification
- This meter features advanced ISE capabilities to work with all ISE testing applications, even the most complex
  - Up to 6 point ISE calibration with timed end point and point-to-point or linear regression slope
  - Non-linear selectable auto-blank and low-concentration range stability
  - Incremental techniques for known or analyte addition and subtraction
  - Adjustable isopoint for temperature-compensated ISE results
- Fix errors without a complete recalibration with the calibration editing function
- Anyone can use this meter with its plain language screen prompts, function keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals or just once with single-shot mode
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods per channel allow for customized procedures to differentiate between tests and/or users
- Two USB and a RS232 port makes for easy data transfer and updating meter software
- Includes data analysis software for easy transferring of stored results or real-time viewing
- Increase your productivity two Orion Star stirrer probe ports allow for simultaneous stirring and control for two separate samples at the same time
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source through the included universal power adapter

- IP54-rated housing handles splashes and cleans up easily
- 3 year meter replacement warranty

### Product Specifications - Orion VERSA STAR Meter

|                                   |                       | CE   TUV   EPA   |
|-----------------------------------|-----------------------|--|
| Maximum Number of Modules         |                       | 4  |
| Methods                           |                       | Up to 10   |
| Datalogging                       | Number of<br>Datalogs | 2000 with time and date stamp                            |
|                                   | Log Function          | Manual, AUTO-READ™, timed interval, single timed reading |
| Communication<br>Ports            |                       | One RS232, two USB ports (one client, one host)          |
| Orion Star Stirrer<br>Probe Ports |                       | 2  |
| Power                             |                       | Universal, AC power adapter,<br>90-260 VAC, 50-60 Hz     |

### Product Specifications - pH/ISE Module (VSTAR-ISE)

| Product Specifications – pH/ISE Module (VSTAR-ISE) |                           |   |  |
|--|---------------------------|---|--|
|  |                           | CE   TUV   EPA  |  |
|  | Range                     | -2.000 to 20.000  |  |
|  | Resolution                | 0.1, 0.01, 0.001  |  |
|  | Relative<br>Accuracy      | ±0.002  |  |
| pH   | Calibration<br>Points     | Up to 6   |  |
|  | Calibration<br>Editing    | Yes   |  |
|  | Adjustable<br>ISO Point   | Yes   |  |
|  | Range                     | ±2000.0 mV  |  |
| mV/Rel mV/   | Resolution                | 0.1   |  |
| ORP  | Relative<br>Accuracy      | ±0.2 mV or ±0.05 % of reading whichever is greater  |  |
|  | E <sub>H</sub> ORP Mode   | Yes   |  |
|  | Range                     | 0 to 19999  |  |
|  | Resolution                | Up to 4 significant digits  |  |
|  | Relative<br>Accuracy      | ±0.2 mV or ±0.05 % of reading whichever is greater  |  |
|  | Units                     | ppm, M, mg/L, %, ppb, none  |  |
|  | Calibration<br>Points     | Up to 6   |  |
| ISE  | Calibration<br>Editing    | Yes   |  |
| 102  | Adjustable<br>ISO Point   | Yes   |  |
|  | Calibration<br>Features   | Timed end point, linear point to point,<br>non-linear selectable auto-blank, low<br>concentration range stability     |  |
|  | Incremental<br>Techniques | Single known addition & subtraction,<br>double known addition & subtraction,<br>single analate addition & subtraction |  |
|  | Advanced<br>Features      | Displays electrode type, temperature compensation   |  |
|  | Range                     | -5 to 105 °C, 23 to 221 °F  |  |
|  | Resolution                | 0.1   |  |
| Temperature  | Relative<br>Accuracy      | ±0.1  |  |
|  | Offset<br>Calibration     | 1 point   |  |
| Innute   | pH/ISE                    | BNC, reference pin  |  |
| Inputs   | ATC                       | 8-pin mini-DIN  |  |



That's not all.
We have more to show you.

### **Meter Ordering Information**

|           | _   |  |  |
|-----------|---|--|--|
| Cat. No.  | Description   |  |  |
| VSTAR40A  | Orion VERSA STAR pH/ISE Benchtop Meter Set - Orion VERSA STAR meter - Orion VERSA STAR pH/ISE/temperature module - Electrode arm with redesigned holder - Universal power adapter   |  |  |
| VSTAR40A2 | Orion VERSA STAR pH/ISE Benchtop Meter Kit  Orion VERSA STAR meter  Orion VERSA STAR pH/ISE/temperature module  8102BNUVWP Orion ROSS Ultra pH electrode  096019 Orion Star stirrer probe  927007MD ATC probe, stainless steel  ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle)  Electrode arm with redesigned holder  Universal power adapter                            |  |  |
| VSTAR40B  | Orion VERSA STAR Benchtop Meter with Two pH/ISE Modules Set - Orion VERSA STAR meter - Two Orion VERSA STAR pH/ISE/temperature modules - Electrode arm with redesigned holder - Universal power adapter   |  |  |
| VSTAR40B2 | Orion VERSA STAR Benchtop Meter with Two pH/ISE Modules Kit  Orion VERSA STAR meter  Two Orion VERSA STAR pH/ISE/temperature modules  8102BNUWP Orion ROSS Ultra pH electrode  906019 Orion Star stirrer probe  927007MD ATC probe, stainless steel  ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle)  Two electrode arms with redesigned holders  Universal power adapter |  |  |



### Accessories

| Accessorie | s   |
|------------|---|
| Cat. No.   | Description   |
| VSTAR-PH   | Orion VERSA STAR pH/temperature module                                    |
| VSTAR-ISE  | Orion VERSA STAR pH/ISE/temperature module                                |
| VSTAR-CND  | Orion VERSA STAR conductivity/temperature module                          |
| VSTAR-RD   | Orion VERSA STAR RDO/DO/temperature module                                |
| VSTAR-LR   | Orion VERSA STAR pH/temperature with LogR technology module               |
| STARA-BEA  | Benchtop electrode arm for Orion Star A-series meters                     |
| STARA-HB   | Freestanding base for use with Orion Star A-series benchtop electrode arm |
| STARA-PWR  | Replacement universal power adapter for Orion<br>VERSA STAR meters        |
| 096019     | Orion Star stirrer probe  |
| 8102BNUWP  | Orion ROSS Ultra pH electrode, refillable, glass body                     |
| 8172BNWP   | Orion ROSS Sure-Flow pH electrode, glass body                             |
| 8156BNUWP  | Orion ROSS Ultra refillable pH electrode, epoxy body                      |
| 8157BNUMD  | Orion ROSS Triode pH/ATC electrode, refillable, epoxy body                |
| 8302BNUMD  | Orion ROSS Triode pH/ATC electrode, refillable, glass body                |
| 927007MD   | Orion ATC probe, stainless steel body                                     |
| 927005MD   | Orion ATC probe, epoxy body   |
| 013005MD   | Orion DuraProbe conductivity cell, 4-cell, K=0.475                        |
| 013016MD   | Orion conductivity cell, flow through, K=0.1                              |
| 083005MD   | Orion dissolved oxygen probe, polarographic                               |
| 086030MD   | Orion auto-stir DO/BOD probe, polarographic                               |
| 087010MD   | Orion RDO optical DO probe, 3 m cable                                     |
| 810199     | Orion ROSS all-in-one pH buffer kit                                       |
| 910104     | Orion pH 4.01 buffer, 475 mL  |
| 910107     | Orion pH 7.00 buffer, 475 mL  |
| 910110     | Orion pH 10.01 buffer, 475 mL   |
| 011008     | Orion conductivity standard, 100 $\mu$ S/cm, 5 x 60 mL bottles            |
| 011007     | Orion conductivity standard, 1413 $\mu$ S/cm, 5 x 60 mL bottles           |
| 011006     | Orion conductivity standard, 12.9 mS/cm, 5 x 60 mL bottles                |
| 1010001    | Orion conductivity calibration resistor kit                               |
| 9512HPBNWP | Orion high-performance ammonia ISE  |
| 9512BNWP   | Orion standard ammonia ISE  |
| 951007     | Ammonia standard, 1000 ppm, 475 mL  |
| 951210     | Low-level ammonia ISA   |
| 951211     | Ammonia ISA, 475 mL   |
| 951209     | Fill solution for Orion high-performance ammonia ISE, 60 mL               |
| 951202     | Fill solution for standard ammonia ISE, 60 mL                             |
| 951213     | Ammonia ISE storage solution  |
| 9609BNWP   | Orion fluoride combination ISE  |
| 040906     | Fluoride standard, 1 ppm with TISAB II, 475 mL                            |
| 040907     | Fluoride standard, 2 ppm with TISAB II, 475 mL                            |
| 040908     | Fluoride standard, 10 ppm with TISAB II, 475 mL                           |
| 940909     | TISAB II, 1 gallon  |
| 8611BNWP   | Orion ROSS sodium combination ISE with sodium standards                   |
|            |   |

ISE \_

ORP

Cond

RDO 📉

DO





### Thermo Scientific Orion DUAL STAR™

Dual Channel pH/ISE Benchtop Meter

Display shows results, electrode status, ATC temperature and channel settings. If a method is selected, the number shows as well

Power key turns on/off the meter and adjusts backlighting

Shortcut keys for direct access to commonly used functions




Channel key toggles display to show both channels simultaneously or just one

Single channel display also shows electrode serial number, sample ID, date and time

Function keys update according to the menu for easy selection

Numeric keypad for easy data entry in submenus



- Dual channel with simultaneously displayed measurements
- Low concentration range reading stability
- Known addition and subtraction incremental methods
- 5 language choices: English, Spanish, French, German and Italian
- Meter prompts and shortcut keys
- Thermo Scientific Orion EZ Startup<sup>™</sup> for quick and easy meter set-up
- Thermo Scientific Orion SMART STABILITY<sup>™</sup> and SMART AVERAGING<sup>™</sup> for accurate and fast results

- E<sub>H</sub> ORP calibration capability
- 10 methods stored per channel

   retains calibrations among
   multiple electrodes
- 1000 point datalogging with time and date stamping
- Calibration, high point and low point alarms
- Calibration graph with slope data for on-screen viewing
- Star stirrer probe controlled directly from the meter
- Autosampler compatibility
- 3 year meter replacement warranty



| 1 Toddet opecifications |  |  |  |
|-------------------------|--|--|--|
|                         |  | CE   TUV   EPA   |  |
| рН                      | Range                                    | -2.000 to 19.999   |  |
|                         | Resolution                               | 0.1, 0.01, 0.001   |  |
|                         | Relative<br>Accuracy                     | ±0.002   |  |
| -                       | Calibration<br>Points                    | 1 to 6   |  |
|                         | Auto-Buffer<br>Recognition               | USA/NIST and EURO/DIN  |  |
|                         | Range                                    | 0 to 19999   |  |
|                         | Resolution                               | 1,2 or 3 significant digits  |  |
|                         | Relative<br>Accuracy                     | ±0.2 mV or 0.05 %, whichever is greater                                      |  |
|                         | Units                                    | ppm, ppb, M, % or unitless   |  |
| ISE                     | Technique<br>Options                     | Single known addition and subtraction  Double known addition and subtraction |  |
|                         | Calibration<br>Points                    | 2 to 6   |  |
|                         | Calibration<br>Options                   | Linear, point to point,<br>non-linear selectable auto-bank                   |  |
|                         | Range                                    | ±1999.9  |  |
| mV/                     | Resolution                               | 0.1  |  |
| Rel mV/<br>ORP          | Relative<br>Accuracy                     | ±0.2 or 0.05 %, whichever is greater   |  |
|                         | E <sub>H</sub> Calibration<br>Capability | Yes  |  |
| Temperature             | Range                                    | -5 to 105 °C   |  |
|                         | Resolution                               | 0.1  |  |
|                         | Relative<br>Accuracy                     | ±0.1 °C  |  |

### Features

| 1 Catal Co  |  |  |
|---|--|--|
| 5 adjustable levels                                       |  |  |
| 20 (10 per channel with calibration data)                 |  |  |
| Yes   |  |  |
| 1000 with time/date stamp                                 |  |  |
| RS232, USB  |  |  |
| IP54  |  |  |
| Universal AC adapter (110V/220V/240V) included with meter |  |  |
| 36 months replacement                                     |  |  |
|   |  |  |

### **Meter Ordering Information**

| Cat. No. | Description  |
|----------|--|
| 2115000  | Orion DUAL STAR pH/ISE meter     Stainless steel ATC probe (927007MD)     Swing arm electrode stand (090043)     RS232 cable (1010053)     Universal power adapter     User guide  |
| 2115101  | <ul> <li>Orion DUAL STAR pH/ISE meter</li> <li>8102BNUWP ROSS Ultra® pH electrode</li> <li>Stainless steel ATC probe (927007MD)</li> <li>RS232 cable (1010053)</li> <li>Swing arm electrode stand (090043)</li> <li>Universal power adapter</li> <li>User guide</li> </ul>   |
| 2115102  | - Orion DUAL STAR pH/ISE meter - 8157BNUMD ROSS Ultra pH/ATC Triode™ refillable electrode - RS232 cable (1010053) - Universal power adapter - User guide   |
| 2115001  | - Orion DUAL STAR pH/ISE meter - 8102BNUWP ROSS Ultra pH electrode - Stainless steel ATC probe (927007MD) - RS232 cable (1010053) - Swing arm electrode stand (090043) - Pint of pH 4.01 buffer (910104) - Pint of pH 7.00 buffer (910107) - Pint of pH 10.01 buffer (910110) - Pint of ROSS® electrode storage solution (810001) - Universal power adapter - User guide |

### Accessories

| Cat. No. | Description                                       |
|----------|---|
| 096019   | Thermo Scientific Orion stirrer probe with paddle |
| 910104   | pH 4.01 buffer, 475 mL                            |
| 910107   | pH 7.00 buffer, 475 mL                            |
| 910110   | pH 10.01 buffer, 475 mL                           |











\_\_\_

 $\square$ 





pH/ISE Benchtop Meter

The Thermo Scientific Orion Star A214 pH/ISE Benchtop Meter is the ideal choice for dedicated ISE analysis. Whether you use ISE infrequently or have multiple samples every day, get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows pH/mV/RmV/ISE results and temperature along with electrode status, time, date, sample ID, user ID and calibration points. For advanced users, features such as stability and averaging choices provide additional options.

# Keep turning the pages for more great finds. ROSS Storage Solution H 4.01 Buffer No Debt STAR ALL Debt STAR

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Up to 5 point ISE calibration with advanced features such as timed end point, linear point-to-point, non-linear selectable auto-blank and low-concentration range stability
- Fix errors without a complete recalibration with the calibration editing function
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Mix samples without a stir plate with direct control of the Orion Star stirrer probe (096019, sold separately)
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- IP54-rated housing handles splashes and is wall-mountable if table space is a problem
- 3 year meter replacement warranty



| •           |                         | CE   TUV   FCC  |
|-------------|-------------------------|---|
|             | Range                   | -2.000 to 20.000  |
| На          | Resolution              | 0.1, 0.01, 0.001  |
|             | Relative<br>Accuracy    | ±0.002  |
| <b>.</b>    | Calibration Points      | Up to 5   |
|             | Calibration<br>Editing  | Yes   |
|             | Range                   | ±2000.0 mV  |
| mV/Rel mV/  | Resolution              | 0.1   |
| ORP         | Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|             | E <sub>H</sub> ORP Mode | Yes   |
|             | Range                   | 0 to 19999  |
|             | Resolution              | Up to 3 significant digits  |
|             | Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|             | Units                   | ppm, M, mg/L, %, ppb, none  |
| ISE         | Calibration<br>Points   | Up to 5   |
|             | Calibration<br>Editing  | Yes   |
|             | Calibration<br>Features | Timed end point, linear point to point,<br>non-linear selectable auto-blank,<br>low concentration range stability |
|             | Range                   | -5 to 105 °C, 23 to 221 °F  |
|             | Resolution              | 0.1   |
| Temperature | Relative<br>Accuracy    | ±0.1  |
|             | Offset<br>Calibration   | 1 point   |
| Datalogging | Number of Datalogs      | 2000 with time and date stamp   |
|             | Log Function            | Manual, ready (include AUTO-READ $^{\text{\tiny TM}}$ ), timed  |
| Inputs      | pH Electrode            | BNC, reference pin  |
|             | ATC Probe               | 8-pin mini-DIN  |
| Output      |                         | RS232, USB  |
|             | AC Adapter              | Included – universal, 100-240 VAC   |
| Power       | Battery<br>Power        | Optional – 4 AAs  |
|             | Battery Life            | 800 hrs   |

### Meter Ordering Information

| Cat. No.  | Description   |  |  |
|-----------|---|--|--|
| STARA2140 | Orion Star A214 pH/ISE Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter  |  |  |
| STARA2145 | Orion Star A214 pH/ISE Benchtop Meter Kit  - 8102BNUWP Orion ROSS Ultra pH electrode  - 927007MD ATC probe, stainless steel  - 096019 Orion Star stirrer probe  - 810199 ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle)  - Electrode arm with redesigned holder  - Universal power adapter   |  |  |
| STARA2146 | Orion Star A214 pH/ISE Benchtop Meter Kit for Ammonia - 8102BNUWP Orion ROSS Ultra pH electrode - 9512HPBNWP Orion high-performance Ammonia ISE - 927007MD ATC probe, stainless steel - 996019 Orion Star stirrer probe - 951007 Orion Ammonia 1000 ppm standard, 475 mL - 951210 Orion Low-level Ammonia ISA - 951213 Ammonia ISE storage solution - Electrode arm with redesigned holder - Universal power adapter  |  |  |
| STARA2147 | Orion Star A214 pH/ISE Benchtop Meter Kit for Fluoride - 8102BNUWP Orion ROSS Ultra pH electrode - 9609BNWP Orion fluoride ISE - 927007MD ATC probe, stainless steel - 096019 Orion Star stirrer probe - 040906 Orion fluoride standard, 1 ppm with TISABII, 475 mL - 040907 Orion fluoride standard, 2 ppm with TISABII, 475 mL - 040908 Orion fluoride standard, 10 ppm with TISABII, 475 mL - 940909 Orion TISABII solution, 1 gallon - Electrode arm with redesigned holder - Universal power adapter |  |  |
| STARA2148 | Orion Star A214 pH/ISE Benchtop Meter Kit for Sodium  - 8102BNUWP Orion ROSS Ultra pH electrode  - 8611BNWP Orion ROSS sodium ISE with sodium standards  - 927007MD ATC probe, stainless steel  - 096019 Orion Star stirrer probe  - Electrode arm with redesigned holder  - Universal power adapter  |  |  |

### Accessories

| Cat. No. | Description                                     |
|----------|---|
| 810199   | Orion ROSS All-in-One pH buffer kit             |
| 951007   | Ammonia standard, 1000 ppm, 475 mL              |
| 951210   | Low-level ammonia ISA                           |
| 951211   | Ammonia ISA, 475 mL                             |
| 040906   | Fluoride standard, 1 ppm with TISAB II, 475 mL  |
| 040907   | Fluoride standard, 2 ppm with TISAB II, 475 mL  |
| 040908   | Fluoride standard, 10 ppm with TISAB II, 475 mL |
| 940909   | TISAB II, 1 gallon                              |

рН 📕

ISE

ORP

Cond

RDO

 $\square$ 





### Thermo Scientific Orion VERSA STAR 50

pH/Conductivity Benchtop Meters

Ideal for checking solution conductivity and follow-up pH testing. The large color display is easy to read and easy on the eyes. Prompts make accuracy and repeatability easy for even the most complicated testing. If you need to add more channels for more or obtain more parameters from the same, simply remove the protective individual channel cover and add the new module.

### **Features and Benefits**

- Versatility is the name of the game
  - Customize the display to show only the information you need
  - Change out modules to create the meter you need at any time
- Up to 6 point pH calibration with automatic recognition for USA/NIST and DIN buffers and electrode identification
- Fix errors without a complete recalibration with the calibration editing function
- Conductivity, TDS, salinity or resistivity capabilities
  - Selectable reading reference temperatures of 5, 10, 15, 20 or 25  $^{\circ}\text{C}$
  - Linear, non-linear, nLFu or USP/EP curve options to get the most accurate conductivity result for your sample and temperature
  - Choice of practical salinity or natural sea water curves for salinity readings
  - Linear or ISO/EN 27888 curves for TDS
- Anyone can use this meter with its plain language screen prompts, function keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals or just once with the single-shot mode
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods per channel allow for customized procedures to differentiate between tests and/or users
- Two USB and a RS232 port makes for easy data transfer and updating meter software
- Includes data analysis software for easy transferring of stored results or real-time viewing
- Control sample mixing directly from the meter without a separate stir plate when using with the Orion Star stirrer probe (cat. no. 096019, sold separately)
- Increase your productivity with two stirrer probe ports, add stirrers and modules to analyze two separate samples at the same time

- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source through the included universal power adapter
- IP54-rated housing handles splashes and cleans up easily
- 3 year meter replacement warranty

### **Product Specifications - Orion VERSA STAR Meter**

|                                   |                       | CE   TUV   EPA   |
|-----------------------------------|-----------------------|--|
| Maximum Number of Modules         |                       | 4  |
| Methods                           |                       | Up to 10   |
| Datalogging                       | Number of<br>Datalogs | 2000 with time and date stamp                            |
| Datalogging                       | Log Function          | Manual, AUTO-READ™, timed interval, single timed reading |
| Communication<br>Ports            |                       | One RS232, two USB ports (one client, one host)          |
| Orion Star Stirrer<br>Probe Ports |                       | 2  |
| Power                             |                       | Universal, AC power adapter,<br>90-260 VAC, 50-60 Hz     |

### Product Specifications - pH Module (VSTAR-PH)

| -           |                         |  |                |
|-------------|-------------------------|--|----------------|
|             |                         |  | CE   TUV   EPA |
|             | Range                   | -2.000 to 20.000                                   |                |
|             | Resolution              | 0.1, 0.01, 0.001                                   |                |
|             | Relative<br>Accuracy    | ±0.002   |                |
| pH          | Calibration<br>Points   | Up to 6  |                |
|             | Calibration<br>Editing  | Yes  |                |
|             | Adjustable<br>ISO Point | Yes  |                |
|             | Range                   | ±2000.0 mV   |                |
| mV/Rel mV/  | Resolution              | 0.1  |                |
| ORP         | Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater |                |
|             | E <sub>H</sub> ORP Mode | Yes  |                |
|             | Range                   | -5 to 105 °C, 23 to 221 °F                         |                |
| Temperature | Resolution              | 0.1  |                |
|             | Relative<br>Accuracy    | ±0.1   |                |
|             | Offset<br>Calibration   | 1 point  |                |
| Innuto      | pH/ISE                  | BNC, reference pin                                 |                |
| Inputs      | ATC                     | 8-pin mini-DIN                                     |                |



### Product Specifications - Conductivity Module (VSTAR-CND) Meter Ordering Information

| Product Specifications – Conductivity Module (VSTAR-CNI |                              |  |
|---|------------------------------|--|
|   |                              | CE   TUV   EPA                                     |
|   | Range                        | 0.001 µS to 3000 mS                                |
|   | Resolution                   | 0.001 µS minimum;<br>4 significant figures minimum |
|   | Relative<br>Accuracy         | 0.5 % reading ±1 digit                             |
|   | Reference<br>Temperature     | 5, 10, 15, 20, 25 °C                               |
| Conductivity  | Temperature<br>Compensation  | Linear (0 to 10.0 %/°C), nLn, nLFu, USP/EP         |
|   | Compatible<br>Cell Constants | 0.001 to 199.9                                     |
|   | Calibration<br>Points        | Up to 6 points                                     |
|   | Calibration<br>Editing       | Yes  |
|   | Range                        | 2 ohm to 100 meg-ohm                               |
| Resistivity   | Resolution                   | 2 ohms-cm  |
| noolouvity  | Relative<br>Accuracy         | 0.5 % reading ±1 digit                             |
|   | Туре                         | Practical salinity or natural sea water            |
| Salinity  | Range                        | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water       |
| Samily  | Resolution                   | 0.01   |
|   | Relative<br>Accuracy         | ±0.1   |
|   | Range                        | 0 to 200 ppt                                       |
|   | Resolution                   | 4 significant digits                               |
| TDS   | Relative<br>Accuracy         | ±0.5 % reading ±1 digit                            |
|   | TDS Factor<br>Range          | Linear 0.01 to 10.00                               |
|   | Range                        | -5 to 105 °C, 23 to 221 °F                         |
|   | Resolution                   | 0.1  |
| Temperature   | Relative<br>Accuracy         | ±0.1   |
|   | Offset<br>Calibration        | 1 point  |
| Inputs  |                              | 8-pin mini-DIN                                     |

|             | •  |
|-------------|--|
| VSTAR50     | Orion VERSA STAR pH/Conductivity Benchtop Meter Set - Orion VERSA STAR meter - Orion VERSA STAR pH/femperature module - Orion VERSA STAR conductivity/temperature module - Electrode arm with redesigned holder - Universal power adapter  |
| VSTAR52     | Orion VERSA STAR pH/Conductivity Benchtop Meter Kit  Orion VERSA STAR meter with pH/temperature module Orion VERSA STAR conductivity/temperature module 8157BNUMD Orion ROSS Ultra Triode pH/ATC probe 013005MD Orion DuraProbe conductivity cell, K=0.475 810199 ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle) 0111007 conductivity standard, 1413 µS, 5 x 60 mL Electrode arm with redesigned holder Universal power adapter |
| Accessories | <b>S</b>   |
| Cat. No.    | Description  |
| VSTAR-PH    | Orion VERSA STAR pH/temperature module   |
| VSTAR-ISE   | Orion VERSA STAR pH/ISE/temperature module   |
| VSTAR-CND   | Orion VERSA STAR conductivity/temperature module   |
| VSTAR-RD    | Orion VERSA STAR RDO/DO/temperature module   |
| VSTAR-LR    | Orion VERSA STAR pH/temperature with LogR  |

Cat. No.

Description

# 8-pin mini-DIN

| Cat. No.  | Description   |  |
|-----------|---|--|
| VSTAR-PH  | Orion VERSA STAR pH/temperature module                                    |  |
| VSTAR-ISE | Orion VERSA STAR pH/ISE/temperature module                                |  |
| VSTAR-CND | Orion VERSA STAR conductivity/temperature module                          |  |
| VSTAR-RD  | Orion VERSA STAR RDO/DO/temperature module                                |  |
| VSTAR-LR  | Orion VERSA STAR pH/temperature with LogR technology module               |  |
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                     |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm |  |
| STARA-PWR | Replacement universal power adapter for Orion<br>VERSA STAR meters        |  |
| 096019    | Orion Star stirrer probe  |  |
| 8102BNUWP | Orion ROSS Ultra pH electrode, refillable, glass body                     |  |
| 8172BNWP  | Orion ROSS Sure-Flow pH electrode, glass body                             |  |
| 8156BNUWP | Orion ROSS Ultra refillable pH electrode, epoxy body                      |  |
| 8157BNUMD | Orion ROSS Triode pH/ATC electrode, refillable, epoxy body                |  |
| 8302BNUMD | Orion ROSS Triode pH/ATC electrode, refillable, glass body                |  |
| 927007MD  | Orion ATC probe, stainless steel body                                     |  |
| 927005MD  | Orion ATC probe, epoxy body   |  |
| 013005MD  | Orion DuraProbe conductivity cell, 4-cell, K=0.475                        |  |
| 013016MD  | Orion conductivity cell, flow through, K=0.1                              |  |
| 810199    | Orion ROSS all-in-one pH buffer kit                                       |  |
| 910104    | Orion pH 4.01 buffer, 475 mL  |  |
| 910107    | Orion pH 7.00 buffer, 475 mL  |  |
| 910110    | Orion pH 10.01 buffer, 475 mL   |  |
| 011008    | Orion conductivity standard, 100 $\mu$ S/cm, 5 x 60 mL bottles            |  |
| 011007    | Orion conductivity standard, 1413 $\mu$ S/cm, 5 x 60 mL bottles           |  |
| 011006    | Orion conductivity standard, 12.9 mS/cm, 5 x 60 mL bottles                |  |
| 1010001   | Orion conductivity calibration resistor kit                               |  |
|           |   |  |

Cond

DDO

 $\square$ 





pH/Conductivity Benchtop Meter

The Thermo Scientific Orion Star A215 pH/Conductivity Benchtop Meter is the ideal choice for any pharmacology lab. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows individual pH/mV/RmV or conductivity/TDS/salinity/resistivity results with temperature, or all results at once. Additional helpful information includes electrode status, time, date, sample ID, user ID and calibration. For advanced users, features such as stability and averaging choices provide additional options.





- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers and calibration editing to fix errors without a complete recalibration
- Selectable cell constant allows for use with 2- or 4-cell conductivity cells
- For accurate conductivity results, selectable reading reference temperatures of 5, 10, 15, 20 or 25 °C with linear, non-linear, nLFu or USP/EP curve options
- Choice of practical salinity or natural sea water curves for salinity readings
- Linear or ISO/EN 27888 curves for TDS
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Mix samples without a stir plate with direct control of the Orion Star stirrer probe (096019, sold separately)
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- IP54-rated housing handles splashes and is wall-mountable if table space is a problem
- 3 year meter replacement warranty



| Product Sp   | ecifications                                     |  |
|--------------|--|--|
|              |  | CE   TUV   FCC                                     |
|              | Range  | -2.000 to 20.000                                   |
|              | Resolution                                       | 0.1, 0.01, 0.001                                   |
| рH           | Relative<br>Accuracy                             | ±0.002   |
|              | Calibration<br>Points                            | Up to 5  |
|              | Calibration<br>Editing                           | Yes  |
|              | Range  | ±2000.0 mV   |
| mV/Rel mV/   | Resolution                                       | 0.1  |
| ORP          | Relative   | ±0.2 mV or ±0.05 % of reading                      |
|              | Accuracy<br>E <sub>H</sub> ORP Mode              | whichever is greater Yes                           |
|              | Range  | 0.001 µS to 3000 mS                                |
|              | Resolution                                       | 0.001 µS minimum;<br>4 significant figures minimum |
|              | Relative<br>Accuracy                             | 0.5 % reading ±1 digit                             |
|              | Reference<br>Temperature                         | 5, 10, 15, 20, 25 °C (default)                     |
| Conductivity | Temperature<br>Compensation                      | Linear (0 to 10.0 %/°C),<br>nLn, nLFu, USP/EP      |
|              | Compatible<br>Cell Constants                     | 0.001 to 199.9                                     |
|              | Calibration<br>Points                            | Up to 5 points                                     |
|              | Calibration<br>Editing                           | Yes  |
|              | Range  | 2 ohm to 100 meg-ohm                               |
| Resistivity  | Resolution                                       | 2 ohms-cm  |
| licolouvity  | Relative<br>Accuracy                             | 0.5 % reading ±1 digit                             |
|              | Туре   | Practical salinity or natural sea water            |
| Calinity     | Range  | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water       |
| Salinity     | Resolution                                       | 0.01   |
|              | Relative<br>Accuracy                             | ±0.1   |
|              | Range  | 1 to 200 ppt                                       |
|              | Resolution                                       | 4 significant digits                               |
| TDS          | Relative<br>Accuracy                             | 0.5 % reading ±1 digit                             |
|              | TDS Factor<br>Range                              | Linear 0.01 to 10.00, default 0.49                 |
|              | Range  | -5 to 105 °C, 23 to 221 °F                         |
|              | Resolution                                       | 0.1  |
| Temperature  | Relative<br>Accuracy                             | ±0.1   |
|              | Offset<br>Calibration                            | 1 point  |
| Datalogging  | Number of Datalogs                               | 2000 with time and date stamp                      |
|              | Log Function                                     | Manual, ready (include AUTO-READ™), timed          |
|              | pH Electrode                                     | BNC, reference pin                                 |
| Inputs       | Conductivity with Temperature Probe or ATC Probe | 8-pin mini-DIN                                     |
| Output       |  | RS232, USB   |
|              | AC Adapter                                       | Included – universal, 100-240 VAC                  |
| Power        | Battery Power                                    | Optional – 4 AAs                                   |
|              | Battery Life                                     | 800 hrs  |
|              |  |  |

### **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA2150 | Orion Star A215 pH/Conductivity Benchtop Multiparameter Meter - Electrode arm with redesigned holder - Universal power adapter  |
| STARA2155 | Orion Star A215 pH/Conductivity Benchtop Multiparameter Meter Kit - 815 7 BNUMD Orion ROSS Ultra Triode pH/ATC electrode - 013005 MD Orion DuraProbe conductivity cell - ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle) - 011007 Orion 1413 µS conductivity standard, 5 x 60 mL - Electrode arm with redesigned holder - Universal power adapter |

### Accessories

| Cat. No.  | Description   |  |
|-----------|---|--|
| 096019    | Orion Star stirrer probe  |  |
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                     |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm |  |
| 8302BNUMD | Orion ROSS Triode 3-in-1 pH/ATC probe, refillable, glass body             |  |
| 8102BNUWP | Orion ROSS Ultra refillable pH electrode, glass body                      |  |
| 927007MD  | Orion ATC probe, stainless steel  |  |
| 810199    | Orion ROSS All-in-One pH buffer kit                                       |  |
| 910104    | Orion pH 4.01 buffer, 475 mL  |  |
| 910107    | Orion pH 7.00 buffer, 475 mL  |  |
| 910110    | Orion pH 10.01 buffer, 475 mL   |  |
| 013005MD  | Orion DuraProbe conductivity cell, 4-cell, K=0.475                        |  |
| 013016MD  | Orion 2-electrode conductivity cell with flow cell, K=0.1                 |  |
| 011008    | Orion conductivity standard, 100 $\mu$ S/cm, 5 $\times$ 60 mL bottles     |  |
| 011007    | Orion conductivity standard, 1413 $\mu$ S/cm, 5 $\times$ 60 mL bottles    |  |
| 011006    | Orion conductivity standard, 12.9 mS/cm, $5 \times 60$ mL bottles         |  |
| 1010001   | Orion conductivity calibration resistor kit                               |  |



рН

ISE \_

ORP

Cond

RDO 🔣

DO





pH/RDO/Dissolved Oxygen Benchtop Meter

The Thermo Scientific Orion Star A216 pH/RDO/Dissolved Oxygen Benchtop Meter is great for checking pH and dissolved oxygen levels at the same time from various wastewater samples. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows individual pH/mV/RmV result, dissolved oxygen in percent saturation or concentration results with temperature, or all results at once. Additional helpful information includes electrode status, time, date, sample ID, user ID and calibration. For advanced users, features such as stability and averaging choices provide additional options.



- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers and calibration editing to fix errors without a complete recalibration
- No need to worry about DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- Accurate DO results are easier with automatic temperature, manual salinity and automatic/manual barometric pressure compensation
- Calibrate DO by using water-saturated air, air-saturated water, a Winkler titration or zero point calibration
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- To mix samples without a stir plate or for faster BOD readings, meter can directly control the 096019 Orion Star stirrer or 086030MD Orion AUTO-STIR™ DO probes (each sold separately)
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- IP54-rated housing handles splashes and is wall-mountable if table space is a problem
- 3 year meter replacement warranty



| Product Sp          | ecifications                                   | •   |
|---------------------|--|---|
|                     |  | CE   TUV   FCC  |
|                     | Range  | -2.000 to 20.000  |
|                     | Resolution                                     | 0.1, 0.01, 0.001  |
| рH                  | Relative<br>Accuracy                           | ±0.002  |
|                     | Calibration<br>Points                          | Up to 5   |
|                     | Calibration<br>Editing                         | Yes   |
|                     | Range  | ±2000.0 mV  |
| mW/Dol mW/          | Resolution                                     | 0.1   |
| mV/Rel mV/<br>ORP   | Relative                                       | ±0.2 mV or ±0.05 % of reading   |
|                     | Accuracy                                       | whichever is greater  |
|                     | E <sub>H</sub> ORP Mode                        | Yes   |
|                     | Concentration                                  | - Polarographic   |
|                     | Range  | 0 to 90 mg/L  |
|                     | Resolution                                     | 0.01, 0.1   |
|                     | Relative<br>Accuracy                           | ±0.2  |
|                     | % Saturation -                                 | Polarographic   |
|                     | Range  | 0 to 600  |
|                     | Resolution                                     | 0.1, 1  |
|                     | Relative<br>Accuracy                           | ±2 %  |
|                     | Concentration                                  | - RDO   |
|                     | Range  | 0 to 50 mg/L  |
|                     | Resolution                                     | 0.01, 0.1   |
| Dissolved<br>Oxygen | Relative<br>Accuracy                           | ±0.1 mg/L up to 8 mg/L;<br>±0.2 mg/L from 8 to 20 mg/L;<br>10 % of reading from 20 to 50 mg/L |
|                     | % Saturation -                                 |   |
|                     | Range  | 0 to 500  |
|                     | Resolution                                     | 0.1, 1  |
|                     | Relative<br>Accuracy                           | ±2 %  |
|                     | Manual<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg  |
|                     | Manual<br>Salinity Factor<br>Correction        | 0 to 45 ppt   |
|                     | Calibration<br>Features                        | Water-saturated air, air saturated water, manual (Winkler) and zero point                     |
|                     | Probe<br>Characteristics                       | Polarographic or RDO  |
|                     | Range  | 0 to 50 °C, 32 to 122 °F  |
|                     | Resolution                                     | 0.1   |
| Temperature         | Relative<br>Accuracy                           | ±0.1  |
|                     | Offset<br>Calibration                          | 1 point   |
| Datalogging         | Number<br>of Datalogs                          | 2000 with time and date stamp   |
|                     | Log Function                                   | Manual, ready (includes AUTO-READ™), timed  |
| Inputs              | pH Electrode<br>DO with                        | BNC, reference pin 9-pin mini-DIN   |
|                     | Temp. Probe                                    |   |
| Output              |  | RS232, USB  |
| Power               | AC Adapter                                     | Optional – universal, 100-240 VAC   |
|                     | Battery Power                                  | Optional – 4 AAs  |
|                     | Battery Life                                   | 800 hrs   |
|                     |  |   |

### **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA2160 | Orion Star A216 pH/Dissolved Oxygen Benchtop<br>Multiparameter Meter<br>- Electrode arm with redesigned holder<br>- Universal power adapter   |
| STARA2165 | Orion Star A216 pH/Dissolved Oxygen Benchtop Multiparameter Meter Kit  - 8157BNUMD Orion ROSS Ultra Triode pH/ATC electrode - 083005MD Orion polarographic DO probe with calibration sleeve - 080513 polarographic DO probe maintenance kit - BOD adapter, funnel and stirrer - 810199 ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle) - Electrode arm with redesigned holder - Universal power adapter |

### Accessories

| Cat. No.  | Description  |  |  |
|-----------|--|--|--|
| 096019    | Orion Star stirrer probe   |  |  |
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                                  |  |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm              |  |  |
| 8302BNUMD | Orion ROSS Triode 3-in-1 pH/ATC probe, refillable, glass body                          |  |  |
| 8102BNUWP | Orion ROSS Ultra refillable pH electrode, glass body                                   |  |  |
| 927007MD  | Orion ATC probe, stainless steel   |  |  |
| 810199    | Orion ROSS All-in-One pH buffer kit  |  |  |
| 910104    | Orion pH 4.01 buffer, 475 mL   |  |  |
| 910107    | Orion pH 7.00 buffer, 475 mL   |  |  |
| 910110    | Orion pH 10.01 buffer, 475 mL  |  |  |
| 083005MD  | Orion polarographic DO probe, 1.5m cable   |  |  |
| 086030MD  | Orion polarographic AUTO-STIR DO/BOD probe, 1.5m cable                                 |  |  |
| 087010MD  | Orion RDO optical DO probe, 3m cable with stainless steel guard and calibration sleeve |  |  |
| 970802    | Orion BOD funnel and stirrer for use with polarographic probes                         |  |  |
| 080360    | Orion BOD bottle adapter for use with polarographic probes                             |  |  |





















### Thermo Scientific Orion VERSA STAR 10 & 80

pH and pH with LogR™ Technology Benchtop Meters

Routine pH measurements are quick and easy with the Orion VERSA STAR 10 and 80 pH benchtop meters. The large color display is easy on to read and easy on the eyes. Prompts make accuracy and repeatability easy for even the most complicated testing. If you need to add more channels for more or obtain more parameters from the same, simply remove the protective individual channel cover and add the new module.

The Orion VERSA STAR 10 pH benchtop meter comes complete with a single channel pH module, offering:

- Up to 6 point pH calibration with automatic recognition for USA/NIST and DIN buffers and electrode identification
- The time-saving option to fix errors without a complete recalibration with the calibration editing function

The Orion VERSA STAR 80 pH with LogR technology benchtop meter includes one LogR module, which fits easily into two channel slots. This technology allows for:

- Having all the features of the VERSA STAR 10 pH benchtop meter
- Preventing wasted efforts and time from bad electrodes by offering the most advanced electrode diagnostics
- Using even less of those expensive samples for testing and still getting accurate pH results since LogR technology allows for temperature-compensated pH readings without a temperature sensor

### **Features and Benefits**

- Versatility is the name of the game
  - Customize the display to show only the information you need
  - Change out modules to create the meter you need at any time
- Anyone can use this meter with its plain language screen prompts, function keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals or just once with single-shot mode
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods per channel allow for customized procedures to differentiate between tests and/or users
- Two USB and a RS232 port makes for easy data transfer and updating meter software
- Includes data analysis software for easy transferring of stored results or real-time viewing

- Control sample mixing directly from the meter without a separate stir plate when using with the Orion Star stirrer probe (cat. no. 096019, sold separately)
- Increase your productivity with two Orion stirrer probe ports, add stirrers and modules to analyze two separate samples at the same time
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source through the included universal power adapter
- IP54-rated housing handles splashes and cleans up easily
- 3 year meter replacement warranty

### Product Specifications - Orion VERSA STAR Meter

|                                   |                       | CE   TUV   EPA   |
|-----------------------------------|-----------------------|--|
| Maximum Number of Modules         |                       | 4  |
| Methods                           |                       | Up to 10   |
| Datalogging                       | Number of<br>Datalogs | 2000 with time and date stamp                            |
|                                   | Log Function          | Manual, AUTO-READ™, timed interval, single timed reading |
| Communication<br>Ports            |                       | One RS232, two USB ports (one client, one host)          |
| Orion Star Stirrer<br>Probe Ports |                       | 2  |
| Power                             |                       | Universal, AC power adapter,<br>90-260 VAC, 50-60 Hz     |





### Product Specifications - pH Module (VSTAR-PH)

|             |                         |  | CE   TUV   EPA |
|-------------|-------------------------|--|----------------|
|             | Range                   | -2.000 to 20.000                                   |                |
|             | Resolution              | 0.1, 0.01, 0.001                                   |                |
|             | Relative<br>Accuracy    | ±0.002   |                |
| pH          | Calibration<br>Points   | Up to 6  |                |
|             | Calibration<br>Editing  | Yes  |                |
|             | Adjustable<br>ISO Point | Yes  |                |
|             | Range                   | ±2000.0 mV   |                |
| mV/Rel mV/  | Resolution              | 0.1  |                |
| ORP         | Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater |                |
|             | E <sub>H</sub> ORP Mode | Yes  |                |
|             | Range                   | -5 to 105 °C, 23 to 221 °F                         |                |
|             | Resolution              | 0.1  |                |
| Temperature | Relative<br>Accuracy    | ±0.1   |                |
|             | Offset<br>Calibration   | 1 point  |                |
| Inputs      | pH/ISE                  | BNC, reference pin                                 |                |
| inputs      | ATC                     | 8-pin mini-DIN                                     |                |

| Product Specifications – pH/LogR Module (VSTAR-LR) |                          |  |                |
|--|--------------------------|--|----------------|
|  |                          |  | CE   TUV   EPA |
|  | Range                    | -2.000 to 20.000                                   |                |
|  | Resolution               | 0.1, 0.01, 0.001                                   |                |
|  | Relative<br>Accuracy     | ±0.002   |                |
| рН   | Calibration<br>Points    | Up to 6  |                |
|  | Calibration<br>Editing   | Yes  |                |
|  | Adjustable<br>ISO Point  | Yes  |                |
|  | Range                    | ±2000.0 mV   |                |
| mV/Rel mV/   | Resolution               | 0.1  |                |
| ORP  | Relative<br>Accuracy     | ±0.2 mV or ±0.05 % of reading whichever is greater |                |
|  | E <sub>H</sub> ORP Mode  | Yes  |                |
|  | Range                    | 0 to 100 °C, 32 to 212 °F                          |                |
|  | Resolution               | 0.1  |                |
|  | Relative<br>Accuracy     | ±0.5   |                |
| LogR   | Resistance<br>Range      | 0 to 6000 MΩ                                       |                |
| Temperature  | Resistance<br>Resolution | 0.1 MΩ up to 1999.9 MΩ, 1 MΩ above 1999.9 MΩ       |                |
|  | Offset<br>Calibration    | Up to 3 points                                     |                |
|  | Calibration<br>Options   | ATC probe, NIST reference                          |                |
|  | Range                    | -5 to 105 °C, 23 to 221 °F                         |                |
| Temperature  | Resolution               | 0.1  |                |
| Measurement<br>without LogR<br>In-Use              | Relative<br>Accuracy     | ±0.1   |                |
|  | Offset<br>Calibration    | 1 point  |                |
| Inputs   | pH                       | BNC, reference pin                                 |                |
| iliputs  | ATC                      | 8-pin mini-DIN                                     |                |

### **Meter Ordering Information**

| Cat. No. | Description  |  |  |
|----------|--|--|--|
| VSTAR10  | Orion VERSA STAR pH Benchtop Meter Set - Orion VERSA STAR meter with pH/temperature module - Electrode arm with redesigned holder - Universal power adapter  |  |  |
| VSTAR12  | Orion VERSA STAR pH Benchtop Meter Kit  Orion VERSA STAR meter with pH/temperature module  8302BNUMD Orion ROSS Triode pH/ATC probe  096019 Orion Star stirrer probe  810199 ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle)  Electrode arm with redesigned holder  Universal power adapter  |  |  |
| VSTAR80  | Orion VERSA STAR pH with LogR Benchtop Meter Set - Orion VERSA STAR meter - Orion VERSA STAR pH/temperature with LogR technology module - 927007MD ATC probe, stainless steel - Electrode arm with redesigned holder - Universal power adapter   |  |  |
| VSTAR82  | Orion VERSA STAR pH Benchtop Meter Kit  Orion VERSA STAR meter  Orion VERSA STAR pH/temperature with LogR technology module  8172BNWP Orion ROSS Sure-Flow pH electrode  927007MD ATC probe, stainless steel  ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle)  810007 ROSS fill solution, 5 x 60 mL  Electrode arm with redesigned holder  Universal power adapter |  |  |

### Accessories

| 40003001103 |  |  |
|-------------|--|--|
| Cat. No.    | Description  |  |
| VSTAR-ISE   | Orion VERSA STAR pH/ISE/temperature module                 |  |
| VSTAR-CND   | Orion VERSA STAR conductivity/temperature module           |  |
| VSTAR-RD    | Orion VERSA STAR RDO/DO/temperature module                 |  |
| 096019      | Orion Star stirrer probe                                   |  |
| 8102BNUWP   | Orion ROSS Ultra pH electrode, refillable, glass body      |  |
| 8172BNWP    | Orion ROSS Sure-Flow pH electrode, glass body              |  |
| 8156BNUWP   | Orion ROSS Ultra refillable pH electrode, epoxy body       |  |
| 8157BNUMD   | Orion ROSS Triode pH/ATC electrode, refillable, epoxy body |  |
| 8302BNUMD   | Orion ROSS Triode pH/ATC electrode, refillable, glass body |  |
| 927007MD    | Orion ATC probe, stainless steel body                      |  |
| 927005MD    | Orion ATC probe, epoxy body                                |  |
| 810199      | Orion ROSS all-in-one pH buffer kit                        |  |
| 910104      | Orion pH 4.01 buffer, 475 mL                               |  |
| 910107      | Orion pH 7.00 buffer, 475 mL                               |  |
| 910110      | Orion pH 10.01 buffer, 475 mL                              |  |

Continue flipping through to discover more!









Cond

RDO 📗

DO





pH Benchtop Meter

The Thermo Scientific Orion Star A211 pH Benchtop Meter is the ideal choice for dedicated pH testing in water, process, research and QA/QC samples. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows pH/mV/RmV and temperature along with electrode status, time, date, sample ID, user ID and calibration points. For advanced users, features such as stability and averaging choices provide additional options.



- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Fix calibration errors without a complete recalibration with calibration editing
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Mix samples without a stir plate with direct control of the Orion Star stirrer probe (096019, sold separately)
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- IP54-rated housing handles splashes and is wall-mountable if table space is a problem
- 3 year meter replacement warranty





|             |                         | CE   TUV   FCC                                     |
|-------------|-------------------------|--|
| рН          | Range                   | -2.000 to 20.000                                   |
|             | Resolution              | 0.1, 0.01, 0.001                                   |
|             | Relative<br>Accuracy    | ±0.002   |
|             | Calibration<br>Points   | Up to 5  |
|             | Calibration<br>Editing  | Yes  |
|             | Range                   | ±2000.0 mV   |
| mV/Rel mV/  | Resolution              | 0.1  |
| ORP         | Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater |
|             | E <sub>H</sub> ORP Mode | Yes  |
|             | Range                   | -5 to 105 °C, 23 to 221 °F                         |
|             | Resolution              | 0.1  |
| Temperature | Relative<br>Accuracy    | ±0.1   |
|             | Offset<br>Calibration   | 1 point  |
| Datalogging | Number of Datalogs      | 2000 with time and date stamp                      |
|             | <b>Log Function</b>     | Manual, ready (includes AUTO-READ™), timed         |
| Inputs      | pH Electrode            | BNC, reference pin                                 |
| iliputs     | ATC Probe               | 8-pin mini-DIN                                     |
| Output      |                         | RS232, USB   |
|             | AC Adapter              | Included - universal, 100-240 VAC                  |
| Power       | Battery<br>Power        | Optional – 4 AAs                                   |
|             | Battery Life            | 800 hrs  |

### **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA2110 | Orion Star A211 pH Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter  |
| STARA2115 | Orion Star A211 pH Benchtop Meter Kit  - 8302BNUMD Orion ROSS Triode 3-in-1 pH/ATC probe  - 810199 ROSS solution kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; cleaning solution; and pH electrode storage bottle)  - Electrode arm with redesigned holder  - Universal power adapter   |
| STARA2116 | Orion Star A211 pH Benchtop Meter Kit with pH/ATC Triode Electrode  - 8157BNUMD Orion ROSS Ultra pH/ATC Triode electrode with epoxy body, refillable  - 810199 All-in-One pH Buffer Kit (475 mL each of pH 4, 7 and 10 buffers; storage solution; and pH electrode storage bottle)  - Electrode arm with redesigned holder  - Universal power adapter |

### Accessories

| Cat. No.  | Description   |  |  |
|-----------|---|--|--|
| 096019    | Orion Star stirrer probe  |  |  |
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                     |  |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm |  |  |
| 8302BNUMD | Orion ROSS Triode 3-in-1 pH/ATC probe, refillable, glass body             |  |  |
| 8102BNUWP | Orion ROSS Ultra refillable pH electrode, glass body                      |  |  |
| 927007MD  | Orion ATC probe, stainless steel  |  |  |
| 810199    | Orion ROSS All-in-One pH buffer kit                                       |  |  |
| 910199    | Orion Standard All-in-One pH buffer kit                                   |  |  |
| 910104    | Orion pH 4.01 buffer, 475 mL  |  |  |
| 910107    | Orion pH 7.00 buffer, 475 mL  |  |  |
| 910110    | Orion pH 10.01 buffer, 475 mL   |  |  |









RDO 📗

DO 📗





pH Benchtop Meter

Thermo Scientific Orion Star A111 pH Benchtop Meters combine simplicity with accuracy for dedicated pH applications. A large LCD displays pH or mV readings along with temperature. Icons provide quick updates on battery life, electrode status and calibration information. Simple button layout and onscreen messages help with calibration and setup menu choices.

- Large, informative screen clearly shows the main information, including pH or mV reading, temperature in °C or °F, meter mode, calibration information and power source
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Up to 3 point pH calibration
- Automatic recognition of USA/NIST and DIN buffers
- Easily recall calibration slope data for procedure checks
- Non-volatile memory holds up to 50 data points
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Plug the meter into almost any power outlet with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- The IP54-rated housing can handle splashes and can be wall-mounted if table space is a problem
- 3 year meter replacement warranty







|             |                                      | CE   TUV   | FCC |
|-------------|--------------------------------------|--|-----|
| рН          | Range                                | -2.00 to 16.00                                     |     |
|             | Resolution                           | 0.1, 0.01  |     |
|             | Relative<br>Accuracy                 | ±0.01  |     |
|             | pH<br>Calibration<br>Points          | Up to 3  |     |
|             | Range - mV                           | ±1600.0 mV   |     |
|             | Range - RmV                          | ±1999.9 mV   |     |
| mV/Rel mV/  | Resolution                           | 0.1  |     |
| ORP         | Relative<br>Accuracy                 | ±0.2 mV or ±0.05 % of reading whichever is greater |     |
|             | E <sub>H</sub> ORP Mode              | No   |     |
|             | Range                                | -5 to 105 °C, 23 to 221 °F                         |     |
|             | Resolution                           | 0.1  |     |
| T           | Relative<br>Accuracy                 | ±0.1   |     |
| Temperature | Temperature<br>Compensation          | Automatic or manual                                |     |
|             | Temperature<br>Offset<br>Calibration | Yes  |     |
| Datalogging | Number of Datalogs                   | 50   |     |
|             | <b>Log Function</b>                  | Manual, Automatic with AUTO-READ™                  |     |
| Inputs      | pH Electrode                         | BNC  |     |
| присо       | ATC Probe                            | 8-pin mini-DIN                                     |     |
|             | AC Adapter                           | Included – universal, 100-240 VAC                  |     |
| Power       | Battery<br>Power                     | Optional – 4 AAs                                   |     |
|             | Battery Life                         | 2000 hrs   |     |

### **Meter Ordering Information**

| Cat. No.  | Description  |
|-----------|--|
| STARA1110 | Orion Star A111 pH Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter   |
| STARA1115 | Orion Star A111 pH Benchtop Meter Kit  - 9157BNMD Orion Triode 3-in-1 pH/ATC probe  - Solution kit (60 mL each of pH 4, 7 and 10 buffers; storage solution and cleaning solution)  - Electrode arm with redesigned holder  - Universal power adapter                     |
| STARA1116 | Orion Star A111 pH Benchtop Meter Kit with TRIS-compatible Electrode  - GD9156BNWP Orion Green combination double-junction refillable pH electrode with epoxy body - 927007MD stainless steel ATC probe - Electrode arm with redesigned holder - Universal power adapter |

### Accessories

| Cat. No.  | Description  |  |  |
|-----------|--|--|--|
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters  |  |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm  |  |  |
| 9157BNMD  | Orion Triode 3-in-1 pH/ATC probe, refillable, epoxy body   |  |  |
| 9107BNMD  | Orion Triode 3-in-1 pH/ATC probe, gel-filled, epoxy body   |  |  |
| 910199    | Orion Standard All-in-One pH buffer kit - 475 mL each of 4.01, 7.00 and 10.01 buffers and pH electrode storage solution with pH electrode storage bottle |  |  |
| 810199    | Orion ROSS All-in-One pH buffer kit  |  |  |
| 910104    | Orion pH 4.01 buffer, 475 mL   |  |  |
| 910107    | Orion pH 7.00 buffer, 475 mL   |  |  |
| 910110    | Orion pH 10.01 buffer, 475 mL  |  |  |
| 910001    | Orion pH electrode storage solution, 475 mL  |  |  |
| 916099    | Orion pH buffer assortment pack – 60 mL each of pH 4, 7 and 10 buffers; storage solution and cleaning solution D   |  |  |



рН 🔣

ISE \_

ORP

Cond

RDO 📗

DO **I** 





### **Thermo Scientific Orion VERSA STAR 20**

Conductivity Benchtop Meters

Conductivity tests can be as clear as ultrapure water samples. The large color display is easy on to read and easy on the eyes. Prompts make accuracy and repeatability easy for even the most complicated testing. If you need to add more channels for more or obtain more parameters from the same, simply remove the protective individual channel cover and add the new module.



- Versatility is the name of the game
  - Customize the display to show only the information you need
- Change modules to create the meter you need at any time
- Selectable reading reference temperatures of 5, 10, 15, 20 or 25 °C
- Linear, non-linear, nLFu or USP/EP curve options to get the most accurate conductivity result for your sample and temperature
- Choice of practical salinity or natural sea water curves for salinity readings
- Linear or ISO/EN 27888 curves for TDS
- Anyone can use this meter with its plain language screen prompts, function keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals or just once with single-shot mode
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods per channel allow for customized procedures to differentiate between tests and/or users
- Two USB and a RS232 port makes for easy data transfer and updating meter software
- Includes data analysis software for easy transferring of stored results or real-time viewing
- Control sample mixing directly from the meter without a separate stir plate when using with the Orion Star stirrer probe (cat. no. 096019, sold separately)
- Increase your productivity with two Orion stirrer probe ports, add stirrers and modules to analyze two separate samples at the same time
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source through the included universal power adapter
- IP54-rated housing handles splashes and cleans up easily
- 3 year meter replacement warranty



### Product Specifications - Orion VERSA STAR Meter

|                                   |                       | CE   TUV   EPA   |
|-----------------------------------|-----------------------|--|
| Maximum Number of Modules         |                       | 4  |
| Methods                           |                       | Up to 10   |
| Datalogging                       | Number of<br>Datalogs | 2000 with time and date stamp                            |
|                                   | Log Function          | Manual, AUTO-READ™, timed interval, single timed reading |
| Communication<br>Ports            |                       | One RS232, two USB ports (one client, one host)          |
| Orion Star Stirrer<br>Probe Ports |                       | 2  |
| Power                             |                       | Universal, AC power adapter,<br>90-260 VAC, 50-60 Hz     |

### Product Specifications - Conductivity Module (VSTAR-CND)

| i roddot op  | oomoa nomo                   | - Conductivity Module (VOTAIT-ON                   |
|--------------|------------------------------|--|
|              |                              | CE   TUV   EPA                                     |
|              | Range                        | 0.001 µS to 3000 mS                                |
|              | Resolution                   | 0.001 µS minimum;<br>4 significant figures minimum |
|              | Relative<br>Accuracy         | 0.5 % reading ±1 digit                             |
|              | Reference<br>Temperature     | 5, 10, 15, 20, 25 °C                               |
| Conductivity | Temperature<br>Compensation  | Linear (0 to 10.0 %/°C), nLn, nLFu, USP/EP         |
|              | Compatible<br>Cell Constants | 0.001 to 199.9                                     |
|              | Calibration<br>Points        | Up to 6 points                                     |
|              | Calibration<br>Editing       | Yes  |
|              | Range                        | 2 ohm to 100 meg-ohm                               |
| Resistivity  | Resolution                   | 2 ohms-cm  |
|              | Relative<br>Accuracy         | 0.5 % reading ±1 digit                             |
|              | Туре                         | Practical salinity or natural sea water            |
| Salinity     | Range                        | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water       |
| Samily       | Resolution                   | 0.01   |
|              | Relative<br>Accuracy         | ±0.1   |
|              | Range                        | 0 to 200 ppt                                       |
|              | Resolution                   | 4 significant digits                               |
| TDS          | Relative<br>Accuracy         | ±0.5 % reading ±1 digit                            |
|              | TDS Factor<br>Range          | Linear 0.01 to 10.00                               |
|              | Range                        | -5 to 105 °C, 23 to 221 °F                         |
| Temperature  | Resolution                   | 0.1  |
|              | Relative<br>Accuracy         | ±0.1   |
|              | Offset<br>Calibration        | 1 point  |
| Input        |                              | 8-pin mini-DIN                                     |
|              |                              |  |

### **Meter Ordering Information**

| Cat. No. | Description  |
|----------|--|
| VSTAR20  | Orion VERSA STAR Conductivity Benchtop Meter Set - Orion VERSA STAR meter - Orion VERSA STAR conductivity/temperature module - Electrode arm with redesigned holder - Universal power adapter  |
| VSTAR22  | Orion VERSA STAR Conductivity Benchtop Meter Kit  Orion VERSA STAR meter  Orion VERSA STAR conductivity/temperature module  O13005MD Orion DuraProbe conductivity cell, K=0.475  O11007 conductivity standard, 1413 µS, 5 x 60 mL  Electrode arm with redesigned holder  Universal power adapter |

### Accessories

| Cat. No.  | Description   |
|-----------|---|
| VSTAR-PH  | Orion VERSA STAR pH/temperature module                                    |
| VSTAR-ISE | Orion VERSA STAR pH/ISE/temperature module                                |
| VSTAR-CND | Orion VERSA STAR conductivity/temperature module                          |
| VSTAR-RD  | Orion VERSA STAR RDO/DO/temperature module                                |
| VSTAR-LR  | Orion VERSA STAR pH/temperature with LogR technology module               |
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                     |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm |
| STARA-PWR | Replacement universal power adapter for Orion<br>VERSA STAR meters        |
| 096019    | Orion Star stirrer probe  |
| 013005MD  | Orion DuraProbe conductivity cell, 4-cell, K=0.475                        |
| 013016MD  | Orion conductivity cell, flow through, K=0.1                              |
| 011008    | Orion conductivity standard, 100 µS/cm, 5 x 60 mL bottles                 |
| 011007    | Orion conductivity standard, 1413 µS/cm, 5 x 60 mL bottles                |
| 011006    | Orion conductivity standard, 12.9 mS/cm, 5 x 60 mL bottles                |
| 1010001   | Orion conductivity calibration resistor kit                               |

















Conductivity Benchtop Meter

The Thermo Scientific Orion Star A212 Conductivity Benchtop Meter is the ideal choice to help meet USP requirements. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows conductivity/TDS/salinity/resistivity and temperature along with time, date, sample ID, user ID and calibration information. For advanced users, features such as stability and averaging choices provide additional options.



- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Selectable cell constant allows for use with 2- or 4-cell conductivity cells
- For accurate results, selectable reading reference temperatures of 5, 10, 15, 20 or 25 °C with linear, non-linear, nLFu or USP/EP curve options
- Choice of practical salinity or natural sea water curves for salinity readings
- Linear or ISO/EN 27888 curves for TDS
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Mix samples without a stir plate with direct control of the Orion Star stirrer probe (096019, sold separately)
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- IP54-rated housing handles splashes and is wall-mountable if table space is a problem
- 3 year meter replacement warranty



| Product Sp   | ecifications                                 |  |
|--------------|--|--|
|              | _  | CE   TUV   FCC   |
|              | Range  | 0.001 µS to 3000 mS  |
|              | Resolution                                   | 0.001 µS minimum;<br>4 significant figures minimum               |
|              | Relative<br>Accuracy                         | 0.5 % reading ±1 digit   |
|              | Reference<br>Temperature                     | 5, 10, 15, 20, 25 °C   |
| Conductivity | Temperature<br>Compensation                  | Linear (0 to 10.0 %/°C), nLn, nLFu, USP/EP                       |
|              | Compatible<br>Cell Constants                 | 0.001 to 199.9   |
|              | Calibration<br>Points                        | Up to 5 points   |
|              | Calibration<br>Editing                       | Yes  |
|              | Range  | 2 ohm to 100 meg-ohm   |
| Resistivity  | Resolution                                   | 2 ohms-cm  |
|              | Relative<br>Accuracy                         | 0.5 % reading ±1 digit   |
|              | Туре   | Practical salinity or natural sea water                          |
| Salinity     | Range  | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water                     |
| Samily       | Resolution                                   | 0.01   |
|              | Relative<br>Accuracy                         | ±0.1   |
|              | Range  | 1 to 200 ppt   |
|              | Resolution                                   | 4 significant digits   |
| TDS          | Relative<br>Accuracy                         | ±0.5 % reading ±1 digit  |
|              | TDS Factor<br>Range                          | Linear, 0.01 to 10.00  |
|              | Range  | -5 to 105 °C, 23 to 221 °F                                       |
|              | Resolution                                   | 0.1  |
| Temperature  | Relative<br>Accuracy                         | ±0.1   |
|              | Offset<br>Calibration                        | 1 point  |
| Datalogging  | Number<br>of Datalogs                        | 2000 with time and date stamp                                    |
|              | Log Function                                 | Manual, ready (includes AUTO-READ $^{\text{\tiny{TM}}}$ ), timed |
| Inputs       | Conductivity<br>with<br>Temperature<br>Probe | 8-pin mini-DIN   |
| Output       |  | RS232, USB   |
|              | AC Adapter                                   | Included – universal, 100-240 VAC                                |
| Power        | <b>Battery Power</b>                         | Optional – 4 AAs   |
|              | Battery Life                                 | 800 hrs  |

### **Meter Ordering Information**

| Cat. No.  | Description  |  |  |
|-----------|--|--|--|
| STARA2120 | Orion Star A212 Conductivity Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter   |  |  |
| STARA2125 | Orion Star A212 Conductivity Benchtop Meter Kit  - 013005MD Orion DuraProbe conductivity probe, K=0.475  - 011007 conductivity standard, 1413 µS, 5 x 60 mL  - Electrode arm with redesigned holder  - Universal power adapter                               |  |  |
| STARA2126 | Orion Star A212 Conductivity Benchtop Meter Kit for Ultrapure Water  - 013016MD Orion conductivity probe with flow-through cell, K=0.1  - 011008 conductivity standard, 100 µS, 5 x 60 mL  - Electrode arm with redesigned holder  - Universal power adapter |  |  |

### Accessories

| Cat. No.  | Description   |  |  |
|-----------|---|--|--|
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                     |  |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm |  |  |
| 013005MD  | Orion DuraProbe conductivity probe, 4-cell, K=0.475                       |  |  |
| 013016MD  | Orion 2-electrode conductivity probe with flow cell, K=0.1                |  |  |
| 011008    | Orion conductivity standard, 100 μS/cm, 5 x 60 mL bottles                 |  |  |
| 011007    | Orion conductivity standard, 1413 µS/cm, 5 x 60 mL bottles                |  |  |
| 011006    | Orion conductivity standard, 12.9 mS/cm, 5 x 60 mL bottles                |  |  |
| 1010001   | Orion conductivity calibration resistor kit                               |  |  |



рН 📗

ISE \_

ORP \_

Cond

RDO 🔣

DO





Conductivity Benchtop Meter

Thermo Scientific Orion Star A112 Conductivity Benchtop Meters combine simplicity with accuracy for dedicated conductivity measurements. A large LCD displays conductivity or TDS readings along with temperature. Icons provide quick updates on battery life and completed calibration status. Simple button layout and onscreen messages help with calibration and setup menu choices.

- Large, informative screen clearly shows the main information, including conductivity or TDS reading, temperature in °C or °F, meter mode and power source
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Selectable cell constant
- Selectable reading reference temperatures of 20 °C or 25 °C for an accurate result
- Included electrode arm and newly-designed holder make it easier to maintain and place probes into samples
- Plug the meter into almost any power outlet with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- Non-volatile memory holds up to 50 data points
- The IP54-rated housing can handle splashes and can be wall-mounted if table space is a problem
- 3 year meter replacement warranty





| Product Specifications |  |   |
|------------------------|--|---|
|                        |  | CE   TUV   FCC                                    |
| Conductivity           | Range  | 0.1 µS/cm to 200 mS/cm                            |
|                        | Resolution                                   | 0.01 µS minimum;<br>3 significant figures minimum |
|                        | Relative<br>Accuracy                         | 0.5 % reading ±1 digit                            |
|                        | Reference<br>Temperature                     | 20, 25 °C (default)                               |
|                        | Temperature<br>Compensation                  | Linear (0.0 to 10.0 %/°C)                         |
|                        | Compatible<br>Cell Constants                 | 0.001 to 10                                       |
|                        | Calibration<br>Points                        | 1   |
|                        | Range  | 1 to 19,999 mg/L                                  |
|                        | Resolution                                   | 4 significant digits                              |
| TDS                    | Relative<br>Accuracy                         | 0.5 % reading ±1 digit                            |
|                        | TDS Factor<br>Range                          | Linear, 0.01 to 10.00, default 0.49               |
|                        | Range  | -5 to 105 °C, 23 to 221 °F                        |
|                        | Resolution                                   | 0.1   |
| Temperature            | Relative<br>Accuracy                         | ±0.1  |
|                        | Offset<br>Calibration                        | Yes   |
| Datalogging            | Number of Datalogs                           | 50  |
|                        | <b>Log Function</b>                          | Manual, Automatic with AUTO-READ™                 |
| Inputs                 | Conductivity<br>with<br>Temperature<br>Probe | 8-pin mini-DIN                                    |
|                        | AC Adapter                                   | Included – universal, 100-240 VAC                 |
| Power                  | <b>Battery Power</b>                         | Optional – 4 AAs                                  |
|                        | Battery Life                                 | 2000 hrs  |
|                        |  |   |

### **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA1120 | Orion Star A112 Conductivity Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter  |
| STARA1125 | Orion Star A112 Conductivity Benchtop Meter Kit - 011050MD orion conductivity cell, K=1.0 - 011007 conductivity standard, 1413 µS, 5 x 60 mL - Electrode arm with redesigned holder - Universal power adapter |

### Accessories

| Cat. No.  | Description   |  |  |
|-----------|---|--|--|
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                     |  |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm |  |  |
| 011050MD  | Orion 2-electrode conductivity cell, K=1.0                                |  |  |
| 013016MD  | Orion 2-electrode conductivity cell with flow cell, K=0.1                 |  |  |
| 011008    | Orion conductivity standard, 100 $\mu$ S/cm, 5 $\times$ 60 mL bottles     |  |  |
| 011007    | Orion conductivity standard, 1413 $\mu$ S/cm, 5 x 60 mL bottles           |  |  |
| 011006    | Orion conductivity standard, 12.9 mS/cm, 5 x 60 mL bottles                |  |  |
| 1010001   | Orion conductivity calibration resistor kit                               |  |  |



рΗ

ISF

ORP

Cond

RDO 📗

DO





### Thermo Scientific Orion VERSA STAR 30

RDO/Dissolved Oxygen Benchtop Meters

BOD testing is complicated enough — no need for the meter to add to the confusion. The large color display is easy on to read and easy on the eyes. Prompts make accuracy and repeatability easy for even the most complicated testing. If you need to add more channels for more or obtain more parameters from the same, simply remove the protective individual channel cover and add the new module.



- Versatility is the name of the game
  - Customize the display to show only the information you need
  - Change modules to create the meter you need at any time
- Takes the worry out of DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- Compensates automatically for temperature, salinity with manual input and barometric pressure (automatic or manual)
- Allows for calibration using water-saturated air, air-saturated water, a Winkler titration or zero point
- Anyone can use this meter with its plain language screen prompts, function keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals or just once with single-shot mode
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods per channel allow for customized procedures to differentiate between tests and/or users
- Two USB and a RS232 port makes for easy data transfer and updating meter software
- Complimentary data analysis software for easy transferring of stored results, real-time viewing and calculation of BOD, DOUR and SOUR results
- Control sample mixing directly from the meter without a separate stir plate when using with the Orion Star stirrer probe (cat. no. 096019, sold separately)
- Increase your productivity with two Orion Star stirrer probe ports, add stirrers and modules to analyze two separate samples at the same time
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source through the included universal power adapter
- IP54-rated housing handles splashes and cleans up easily
- 3 year meter replacement warranty



### Product Specifications - Orion VERSA STAR Meter

|                                   |                       | CE   TUV   EPA   |
|-----------------------------------|-----------------------|--|
| Maximum Number of Modules         |                       | 4  |
| Methods                           |                       | Up to 10   |
| Datalogging                       | Number of<br>Datalogs | 2000 with time and date stamp                            |
|                                   | Log Function          | Manual, AUTO-READ™, timed interval, single timed reading |
| Communication<br>Ports            |                       | One RS232, two USB ports (one client, one host)          |
| Orion Star Stirrer<br>Probe Ports |                       | 2  |
| Power                             |                       | Universal, AC power adapter,<br>90-260 VAC, 50-60 Hz     |

### Product Specifications - RDO/DO Module (VSTAR-RD)

| Product Specifications – RDO/DO Module (VSTAR-RD) |   |   |
|---|---|---|
|   |   | CE   TUV   EPA  |
|   | Concentration                                     | - Polarographic   |
|   | Range   | 0 to 90 mg/L  |
|   | Resolution  | 0.01, 0.1   |
|   | Relative<br>Accuracy                              | ±0.2  |
|   | % Saturation -                                    | Polarographic   |
|   | Range   | 0 to 600  |
|   | Resolution  | 0.1, 1  |
|   | Relative<br>Accuracy                              | ±2 %  |
|   | Concentration                                     | - RD0   |
|   | Range   | 0 to 50 mg/L  |
|   | Resolution  | 0.01, 0.1   |
| Dissolved<br>Oxygen                               | Relative<br>Accuracy                              | ±0.1 mg/L up to 8 mg/L;<br>±0.2 mg/L from 8 to 20 mg/L;<br>10 % of reading, 20 to 50 mg/L |
|   | % Saturation -                                    | RDO   |
|   | Range   | 0 to 500  |
|   | Resolution  | 0.1, 1  |
|   | Relative<br>Accuracy                              | ±2 %  |
|   | Automatic<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg  |
|   | Salinity<br>Factor<br>Correction                  | Automatic when used with conductivity module and probe or manual; 0 to 45 ppt             |
|   | Calibration<br>Features                           | Water saturated air, air saturated water, manual (Winkler) and zero point                 |
|   | Probe Type  | Polarographic or RDO  |
|   | Range   | 0 to 50 °C, 32 to 122 °F  |
|   | Resolution  | 0.1   |
| Temperature                                       | Relative<br>Accuracy                              | ±0.1  |
|   | Offset<br>Calibration                             | 1 point   |
| Input   |   | 9-pin mini-DIN  |
|   |   |   |

### **Meter Ordering Information**

| Cat. No. | Description   |
|----------|---|
| VSTAR30  | Orion VERSA STAR RDO/Dissolved Oxygen Benchtop Meter Set - Orion VERSA STAR meter - Orion VERSA STAR dissolved oxygen/ temperature module - Electrode arm with redesigned holder - Universal power adapter  |
| VSTAR32  | Orion VERSA STAR RDO/Dissolved Oxygen Benchtop Meter Kit Orion VERSA STAR meter Orion VERSA STAR dissolved oxygen/ temperature module 083005MD Orion polarographic DO probe 080017 calibration sleeve for DO probe 080513 DO probe maintenance kit BOD adapter, funnel and stirrer Electrode arm with redesigned holder Universal power adapter |

### Accessories

| Cat. No.  | Description   |
|-----------|---|
| VSTAR-PH  | Orion VERSA STAR pH/temperature module  |
| VSTAR-ISE | Orion VERSA STAR pH/ISE/temperature module  |
| VSTAR-CND | Orion VERSA STAR conductivity/temperature module  |
| VSTAR-RD  | Orion VERSA STAR RDO/DO/temperature module  |
| VSTAR-LR  | Orion VERSA STAR pH/temperature with LogR technology module   |
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters   |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm   |
| STARA-PWR | Replacement universal power adapter for Orion<br>VERSA STAR meters  |
| 096019    | Orion Star stirrer probe  |
| 083005MD  | Orion dissolved oxygen probe, polarographic   |
| 086030MD  | Orion auto-stir DO/BOD probe, polarographic   |
| 087010MD  | Orion RDO optical DO probe, 3 m cable   |
| 080513    | Maintenance kit, includes electrolyte solution,<br>polishing disk, and two membrane caps for 083005MD<br>and 086030MD |
| 080514    | Polarographic electrolyte solution for 083005MD and 086030MD  |
| 086021A   | Replacement calibration chamber/stand for 086030MD  |
| 087001    | Replacement optical/luminescence-based cap  |
|           |   |

рН 📉

ISE \_

ORP \_

Cond

RDO 🗾

DO 🔣





RDO/Dissolved Oxygen Benchtop Meter

The Thermo Scientific Orion Star A213 RDO/Dissolved Oxygen (DO) Benchtop Meter is the ideal choice to ensure proper oxygen levels. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows dissolved oxygen in percent saturation or concentration and temperature along with time, date, sample ID, user ID and calibration information. For advanced users, features such as stability and averaging choices provide additional options.



- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- No need to worry about DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- Accurate results are easy with automatic temperature, manual salinity and your choice of automatic or manual barometric pressure compensation
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Calibrate using water-saturated air, air-saturated water, a Winkler titration or zero point calibration
- Non-volatile memory holds up to 2000 data points with time and date stamp
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- For faster BOD readings, control the Orion AUTO-STIR<sup>™</sup> polarographic DO probe (086030MD, sold separately) directly from the meter
- Included electrode arm and newly-designed probe holder make it easier to maintain and place probes into samples
- Works with almost every AC power source with the included universal power adapter, or use four AA batteries (sold separately) to run the meter on DC power
- IP54-rated housing handles splashes and is wall-mountable if table space is a problem
- 3 year meter replacement warranty



| r roduct op         | ecilications                                      |   |
|---------------------|---|---|
|                     |   | CE   TUV   FCC  |
|                     |   | - Polarographic   |
|                     | Range   | 0 to 90 mg/L  |
|                     | Resolution  | 0.01, 0.1   |
|                     | Relative<br>Accuracy                              | ±0.2  |
|                     | % Saturation -                                    | Polarographic   |
|                     | Range   | 0 to 600  |
|                     | Resolution  | 0.1, 1  |
|                     | Relative<br>Accuracy                              | ±2 %  |
|                     | Concentration                                     | - RDO   |
|                     | Range   | 0 to 50 mg/L  |
|                     | Resolution  | 0.01, 0.1   |
| Dissolved<br>Oxygen | Relative<br>Accuracy                              | ±0.1 mg/L up to 8 mg/L;<br>±0.2 mg/L from 8 to 20 mg/L;<br>10 % of reading from 20 to 50 mg/L |
| ,,,,                | % Saturation -                                    | RDO   |
|                     | Range   | 0 to 500  |
|                     | Resolution  | 0.1, 1  |
|                     | Relative<br>Accuracy                              | ±2 %  |
|                     | Automatic<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg  |
|                     | Manual<br>Salinity Factor<br>Correction           | 0 to 45 ppt   |
|                     | Calibration<br>Features                           | Water-saturated air, air saturated water,<br>manual (Winkler) and zero point                  |
|                     | Probe<br>Characteristics                          | Polarographic or RDO  |
|                     | Range   | 0 to 50 °C, 32 to 122 °F  |
|                     | Resolution  | 0.1   |
| Temperature         | Relative<br>Accuracy                              | ±0.1  |
|                     | Offset<br>Calibration                             | 1 point   |
| Datalogging         | Number of Datalogs                                | 2000 with time and date stamp   |
|                     | Log Function                                      | Manual, ready (includes AUTO-READ™), timed  |
| Inputs              | DO with<br>Temperature<br>Probe                   | 9-pin mini-DIN  |
| Output              |   | RS232, USB  |
|                     | AC Adapter  | Included – universal, 100-240 VAC   |
| Power               | Battery Power                                     | Optional – 4 AAs  |
|                     | Battery Life                                      | 800 hrs   |
|                     | •   |   |

### **Meter Ordering Information**

| Cat. No.  | Description   |  |  |
|-----------|---|--|--|
| STARA2130 | Orion Star A213 Dissolved Oxygen Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter  |  |  |
| STARA2135 | Orion Star A213 Dissolved Oxygen Benchtop Meter Kit  - 083005MD Orion polarographic DO probe with calibration sleeve  - 080513 polarographic DO probe maintenance kit  - BOD adapter, funnel and stirrer  - Electrode arm with redesigned holder  - Universal power adapter                     |  |  |
| STARA2136 | Orion Star A213 Dissolved Oxygen Benchtop Meter Kit with Self-Stirring Probe  - 086030MD Orion polarographic AUTO-STIR DO/BOD probe  - 080513 DO probe maintenance kit  - 080514 polarographic DO probe electrolyte solution  - Electrode arm with redesigned holder  - Universal power adapter |  |  |

### Accessories

| Cat. No.  | Description  |  |  |
|-----------|--|--|--|
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                                  |  |  |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm              |  |  |
| 083005MD  | Orion polarographic DO probe, 1.5m cable   |  |  |
| 086030MD  | Orion polarographic AUTO-STIR DO/BOD probe,<br>1.5m cable                              |  |  |
| 087010MD  | Orion RDO optical DO probe, 3m cable with stainless steel guard and calibration sleeve |  |  |
| 080513    | Orion polarographic DO probe maintenance kit   |  |  |
| 080514    | Orion polarographic DO probe electrolyte solution, 60 mL                               |  |  |
| 080017    | Calibration sleeve for Orion polarographic DO probes                                   |  |  |
| 087003    | Replacement calibration sleeve for Orion RDO optical probes                            |  |  |
| 087001    | Replacement optical cap for Orion RDO probes   |  |  |
| 970802    | Orion BOD funnel and stirrer for use with polarographic probes                         |  |  |
| 080360    | Orion BOD bottle adapter for use with polarographic probes                             |  |  |







Cond

RDO 🗾

DO 🗾





Dissolved Oxygen Benchtop Meter

Thermo Scientific Orion Star A113 Dissolved Oxygen Benchtop Meters combine simplicity with accuracy for dedicated DO testing. Icons provide quick updates on battery life and completed calibration status. Simple button layout and onscreen messages help with calibration and setup menu choices.

- Large, informative screen clearly shows the main information, including DO reading in mg/L or % saturation, temperature in °C or °F, meter mode and power source
- Automatic temperature, manual salinity input and manual barometric pressure compensation for an accurate result
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Two calibration options to meet your requirements: water-saturated air method or a Winkler titration
- Accepts Thermo Scientific Orion polarographic dissolved oxygen probes
- Non-volatile memory holds up to 50 data points
- Included electrode arm and newly-designed holder make it easier to maintain and place probes into samples
- Plug the meter into almost any power outlet with the included universal power adapter, or use four AA batteries (sold separately)
- The IP54-rated housing can handle splashes and can be wall-mounted if table space is a problem
- 3 year meter replacement warranty





| Product Specifications |  |                                       |  |
|------------------------|--|---------------------------------------|--|
|                        |  | CE   TUV   FCC                        |  |
|                        | Concentration                                  |                                       |  |
|                        | Range  | 0 to 20 mg/L                          |  |
|                        | Resolution                                     | 0.01, 0.1                             |  |
|                        | Relative<br>Accuracy                           | ±0.2                                  |  |
|                        | % Saturation                                   |                                       |  |
|                        | Range  | 0 to 200                              |  |
|                        | Resolution                                     | 0.1, 1                                |  |
| Dissolved<br>Oxygen    | Relative<br>Accuracy                           | ±2 %                                  |  |
|                        | Manual<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg                  |  |
|                        | Manual<br>Salinity Factor<br>Correction        | 0 to 45 ppt                           |  |
|                        | Calibration<br>Features                        | Water-saturated air, manual (Winkler) |  |
|                        | Sensor Type                                    | Polarographic                         |  |
|                        | Range  | 0 to 50 °C, 32 to 122 °F              |  |
|                        | Resolution                                     | 0.1                                   |  |
| Temperature            | Relative<br>Accuracy                           | ±0.1                                  |  |
|                        | Offset<br>Calibration                          | Yes                                   |  |
| Datalogging            | Number of Datalogs                             | 50                                    |  |
|                        | Log Function                                   | Manual, Automatic with AUTO-READ™     |  |
| Inputs                 | DO with<br>Temperature<br>Probe                | 9-pin mini-DIN                        |  |
|                        | Power Cord                                     | Included – universal, 100-240 VAC     |  |
| Power                  | <b>Battery Power</b>                           | Optional – 4 AAs                      |  |
|                        | Battery Life                                   | 2000 hrs                              |  |
|                        |  |                                       |  |

### **Meter Ordering Information**

| Cat. No.  | Description  |
|-----------|--|
| STARA1130 | Orion Star A113 Dissolved Oxygen Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter   |
| STARA1135 | Orion Star A113 Dissolved Oxygen Benchtop Meter Kit  - 083005MD Orion polarographic DO probe  - 080017 calibration sleeve for DO probe  - 080513 DO probe maintenance kit  - Electrode arm with redesigned holder  - Universal power adapter |

### Accessories

| Cat. No.  | Description   |
|-----------|---|
| STARA-BEA | Benchtop electrode arm for Orion Star A-series meters                     |
| STARA-HB  | Freestanding base for use with Orion Star A-series benchtop electrode arm |
| 083005MD  | Orion polarographic DO probe, 1.5m cable                                  |
| 080513    | Orion polarographic DO probe maintenance kit                              |
| 080017    | Calibration sleeve for Orion polarographic DO probes                      |
| 970802    | Orion BOD funnel and stirrer for use with polarographic probes            |
| 080360    | Orion BOD bottle adapter for use with polarographic probes                |

For accurate testing outside the lab, check out the Orion Star A portable meters on the following pages.

рН

ISE

ORP |

Cond

RDO

DO 📉





### Thermo Scientific Orion VERSA STAR Advanced Electrochemistry Benchtop Meter & Modules

Don't see the benchtop meter parameter combination you need for your application? No problem — simply order the Orion VERSA STAR meter and choose the modules you need.

With versatility being the name of the game, you can also choose modules to add or change your existing Orion VERSA STAR meter to meet your latest testing needs. Simply swap out the modules when you need to and no need to wait to power down. The meter automatically recognizes the changes and datalogs remain unchanged and stored in the meter.

Each module includes the necessary probe connections for the measurement parameter. The meter includes two Orion Star stirrer probe ports, one RS232 port, two USB ports, power port and grounding port. Included with the meter are the benchtop electrode arm and individual channel covers to protect the channels not in use.

Thermo

### Meter and Module Ordering Information

| Cat. No.  | Description  |
|-----------|--|
| VSTAR00   | Orion VERSA STAR Benchtop Meter - Electrode arm with redesigned holder - Universal power adapter |
| VSTAR-PH  | Orion VERSA STAR pH/temperature module   |
| VSTAR-ISE | Orion VERSA STAR pH/ISE/temperature module   |
| VSTAR-CND | Orion VERSA STAR conductivity/temperature module   |
| VSTAR-RD  | Orion VERSA STAR RDO/DO/temperature module   |
| VSTAR-LR  | Orion VERSA STAR pH/temperature with LogR technology module                                      |
| STARA-BEA | Benchtop electrode arm   |
| STARA-HB  | Free-standing, weighted base for use with benchtop electrode arm (sold separately)               |



ORION VERSASTAR



ISE

ORP

Cond

RDO

DO







# Thermo Scientific Orion Star A320, A220 and A120-Series Portable Meters

Outside the laboratory, testing brings enough of its own challenges from getting the sample to taking the reading. Just because you test outside the laboratory doesn't mean your testing accuracy or quality needs to suffer. The new Orion Star A series of portable meters features a variety of advantages to handle testing in difficult environments.



Do you have access to an AC power source while testing? Save battery life by running the meter off AC power with the optional adapter







pH/ISE/Conductivity/RDO/Dissolved Oxygen Portable Meter

The Thermo Scientific Orion Star A329 pH/ISE/Conductivity/RDO/Dissolved Oxygen Portable Meter is the top choice for advanced multiparameter measurements on location. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display can show each parameter with temperature individually or all three with temperature on one screen. Additional helpful icons include electrode status, time, date, sample ID, user ID and calibration points. Advanced features such as stability and averaging options provide additional options allowing for laboratory accuracy everywhere you need it.

## **Features and Benefits**

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Fix pH and ISE calibration errors without a complete recalibration with calibration editing
- No need to worry about DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- It's easy to get accurate DO results with automatic temperature, automatic/manual barometric pressure and automatic salinity compensation when used with simultaneously with conductivity
- Calibrate for DO using water-saturated air, air-saturated water, a Winkler titration or zero point calibration
- Selectable cell constant allows for use with 2- or 4-cell conductivity cells
- For advanced accuracy in conductivity readings, selectable reading reference temperatures of 5, 10, 15, 20 or 25 °C with linear, non-linear, nLFu or USP/EP curve options
- Choice of practical salinity or natural sea water curves for salinity readings
- Linear or ISO/EN 27888 curves for TDS
- Whether you use ISE infrequently or have multiple samples every day — up to 5 point ISE calibration with advanced features such as timed end point, linear point-to-point, non-linear selectable auto-blank and low-concentration range stability
- Non-volatile memory holds up to 5000 data points with time and date stamp

- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- User and sample ID fields for additional data accuracy
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty

#### **Product Specifications**

|              | ecilications                 | CE   TUV   FCC  |
|--------------|------------------------------|---|
|              | Range                        | -2.000 to 20.000  |
| рН           | Resolution                   | 0.1, 0.01, 0.001  |
|              | Relative<br>Accuracy         | ±0.002  |
|              | Calibration<br>Points        | Up to 5   |
|              | Calibration<br>Editing       | Yes   |
|              | Range                        | ±2000.0 mV  |
| mV/Rel mV/   | Resolution                   | 0.1   |
| ORP          | Relative<br>Accuracy         | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|              | E <sub>H</sub> ORP Mode      | Yes   |
|              | Range                        | 0 to 19999  |
|              | Resolution                   | Up to 3 significant digits  |
|              | Relative<br>Accuracy         | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|              | Units                        | ppm, M, mg/L, %, ppb, none  |
| ISE          | Calibration<br>Points        | Up to 5   |
|              | Calibration<br>Editing       | Yes   |
|              | Calibration<br>Features      | Timed end point, linear point to point,<br>non-linear selectable auto-blank, low<br>concentration range stability |
|              | Range                        | 0.001 µS to 3000 mS   |
|              | Resolution                   | 0.001 µS minimum;<br>4 significant figures minimum  |
|              | Relative<br>Accuracy         | 0.5 % reading ±1 digit  |
| Conductivity | Reference<br>Temperature     | 5, 10, 15, 20, 25 °C (default)  |
| Conductivity | Temperature<br>Compensation  | Linear (0 to 10.0 %/°C),<br>nLn, nLFu, USP/EP   |
|              | Compatible<br>Cell Constants | 0.001 to 199.9  |
|              | Calibration<br>Points        | Up to 5 points  |
|              | Calibration<br>Editing       | Yes   |
|              | Range                        | 2 ohm to 100 meg-ohm  |
| Resistivity  | Resolution                   | 2 ohms-cm   |
|              | Relative<br>Accuracy         | 0.5 % reading ±1 digit  |





Please turn the page... this is just the beginning of the portable meters.

## **Product Specifications (Continued)**

| Product Sp          | ecifications                                      | (Continued)  |
|---------------------|---|--|
|                     |   | CE   TUV   FCC   |
| Salinity            | Туре  | Practical salinity or natural sea water                                      |
|                     | Range   | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water                                 |
|                     | Resolution  | 0.01   |
|                     | Relative<br>Accuracy                              | ±0.1   |
|                     | Range   | 1 to 200 ppt   |
|                     | Resolution  | 4 significant digits   |
| TDS                 | Relative<br>Accuracy                              | ±0.5 % reading ±1 digit  |
|                     | TDS Factor<br>Range                               | Linear 0.01 to 10.00, default 0.49   |
|                     | Concentration                                     | - Polarographic  |
|                     | Range   | 0 to 90 mg/L   |
|                     | Resolution  | 0.01, 0.1  |
|                     | Relative<br>Accuracy                              | ±0.2   |
|                     | % Saturation -                                    | Polarographic  |
|                     | Range   | 0 to 600   |
|                     | Resolution  | 0.1, 1   |
|                     | Relative<br>Accuracy                              | ±2 %   |
|                     | Concentration                                     | - RD0  |
|                     | Range   | 0 to 50 mg/L   |
|                     | Resolution  | 0.01, 0.1  |
| Dissolved<br>Oxygen | Relative<br>Accuracy                              | ±0.1 mg/L up to 8 mg/L;<br>±0.2 mg/L from 8 to 20 mg/L;                      |
| Охуусп              | % Saturation -                                    | 10 % of reading from 20 to 50 mg/L   |
|                     | Range   | 0 to 500   |
|                     | Resolution  | 0.1, 1   |
|                     | Relative<br>Accuracy                              | ±2 %   |
|                     | Automatic<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg   |
|                     | Salinity Factor<br>Correction                     | Automatic with separate conductivity probe or manual; 0 to 45 ppt            |
|                     | Calibration<br>Features                           | Water-saturated air, air saturated water,<br>manual (Winkler) and zero point |
|                     | Probe<br>Characteristics                          | Polarographic or RDO   |
|                     | Range   | -5 to 105 °C, 23 to 221 °F<br>(0 to 50 °C, 32 to 122 °F with RDO)            |
|                     | Resolution  | 0.1  |
| Temperature         | Relative<br>Accuracy                              | ±0.1   |
|                     | Offset<br>Calibration                             | 1 point  |
| Datalogging         | Number of Datalogs                                | 5000 with time and date stamp  |
|                     | <b>Log Function</b>                               | Manual, ready (includes AUTO-READ $^{\text{\tiny{TM}}}$ ), timed             |
|                     | pH Electrode<br>Conductivity                      | BNC, reference pin   |
| Inputs              | with Temperature Probe or ATC Probe               | 8-pin mini-DIN   |
|                     | DO with<br>Temperature<br>Probe                   | 9-pin mini-DIN   |
| Output              | 11000   | RS232, USB   |
| Jacput              | AC Adapter  | Optional – universal, 100-240 VAC  |
| Power               | Battery Power                                     | Included – 4 AAs   |
|                     | Battery Life                                      | 800 hrs  |

## **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA3290 | Orion Star A329 pH/ISE/Conductivity/Dissolved Oxygen Portable Multiparameter Meter - Four AA batteries (installed)  |
| STARA3295 | Orion Star A329 pH/ISE/Conductivity/Dissolved Oxygen Portable Multiparameter Meter Kit  - 8107 UWMIMD Orion ROSS Triode 3-in-1 pH/ATC probe, 3m cable  - 013010MD Orion DuraProbe conductivity cell, 3m cable  - 087010MD Orion RDO optical/luminescence-based probe with 3m cable, stainless steel guard, calibration sleeve and optical/luminescence-based cap  - 910410 Orion pH 4.01 buffer, 10 pouches  - 910710 Orion pH 7.00 buffer, 10 pouches  - 911010 Orion pH 10.01 buffer, 10 pouches  - 810001 Orion ROSS electrode storage solution, 475 mL - 01100710 Orion 1413 µS conductivity standard, 10 pouches  - 9111110 rinse solution, 10 pouches  - Protective armor with probe holders  - Hard carrying case  - Four AA batteries (installed) |

#### Accessories

| Cat. No.   | Description  |  |  |
|------------|--|--|--|
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |  |  |
| 1210005    | Soft carrying case for Orion Star, Star Plus and Star<br>A-series portable meters                            |  |  |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |  |  |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |  |  |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star A-series armor   |  |  |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |  |  |
| 8107UWMMD  | Orion ROSS Triode 3-in-1 pH/ATC probe, gel-filled, 3m cable  |  |  |
| 910410     | Orion pH 4.01 buffer pouches, 10/pk  |  |  |
| 910710     | Orion pH 7.00 buffer pouches, 10/pk  |  |  |
| 911010     | Orion pH 10.01 buffer pouches, 10/pk   |  |  |
| 911110     | Orion pH rinse solution pouches, 10/pk   |  |  |
| 013010MD   | Orion DuraProbe conductivity cell with 3m cable, 4-cell, K=0.475   |  |  |
| 013020MD   | Orion DuraProbe conductivity cell with 6m cable, 4-cell, K=0.475   |  |  |
| 013025MD   | Orion DuraProbe conductivity cell with 10m cable, 4-cell, K=0.475  |  |  |
| 011008     | Orion conductivity standard, 100 $\mu$ S/cm, 5 x 60 mL bottles   |  |  |
| 01100710   | Orion conductivity standard pouches, 1413 µS/cm, 10/pk   |  |  |
| 01100610   | Orion conductivity standard pouches, 12.9 mS/cm, 10/pk   |  |  |
| 1010001    | Orion conductivity calibration resistor kit  |  |  |
| 087010MD   | Orion RDO optical DO probe, 3m cable with stainless steel guard and calibration sleeve                       |  |  |
| 087020MD   | Orion RDO optical DO probe, 6m cable with stainless steel guard and calibration sleeve                       |  |  |
| 087030MD   | Orion RDO optical DO probe, 10m cable with stainless steel guard and calibration sleeve                      |  |  |
| 087050MD   | Orion RDO optical DO probe, 15m cable with stainless steel guard and calibration sleeve                      |  |  |
| 087100MD   | Orion RDO optical DO probe, 30m cable with stainless steel guard and calibration sleeve                      |  |  |
| 083005MD   | Orion polarographic DO probe, 1.5m cable   |  |  |
| 080513     | Orion polarographic DO probe maintenance kit   |  |  |
| 080514     | Orion polarographic DO probe electrolyte solution, 60 mL   |  |  |
| 080017     | Calibration sleeve for Orion polarographic DO probes   |  |  |
| 087003     | Replacement calibration sleeve for Orion RDO optical probes  |  |  |
| 087001     | Replacement optical cap for Orion RDO probes   |  |  |

рН 🔃

ISE \_

ORP

Cond

RDO 🗾

DO **I** 





pH/ISE Portable Meter

The Thermo Scientific Orion Star A324 pH/ISE Portable Meters is the top choice for advanced pH/ISE measurements on location. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows pH/mV/RmV and temperature along with electrode status, time, date, sample ID, user ID and calibration points. Advanced features such as stability and averaging options provide additional options allowing for laboratory accuracy everywhere you need it.



- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Whether you use ISE infrequently or have multiple samples every day — up to 5 point ISE calibration with advanced features such as timed end point, linear point-to-point, non-linear selectable auto-blank and low-concentration range stability
- Fix errors without a complete recalibration with the calibration editing function
- Non-volatile memory holds up to 5000 data points with time and date stamp
- User and sample ID fields for additional data accuracy
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty



| Product Sp  | ecincations             |   |
|-------------|-------------------------|---|
|             |                         | CE   TUV   FCC  |
| рН          | Range                   | -2.000 to 20.000  |
|             | Resolution              | 0.1, 0.01, 0.001  |
|             | Relative<br>Accuracy    | ±0.002  |
|             | Calibration<br>Points   | Up to 5   |
|             | Calibration<br>Editing  | Yes   |
|             | Range                   | ±2000.0 mV  |
| mV/Rel mV/  | Resolution              | 0.1   |
| ORP         | Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|             | E <sub>H</sub> ORP Mode | Yes   |
|             | Range                   | 0 to 19999  |
|             | Resolution              | Up to 3 significant digits  |
|             | Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater  |
|             | Units                   | ppm, M, mg/L, %, ppb, none  |
| ISE         | Calibration<br>Points   | Up to 5   |
|             | Calibration<br>Editing  | Yes   |
|             | Calibration<br>Features | Timed end point, linear point to point,<br>non-linear selectable auto-blank,<br>low concentration range stability |
|             | Range                   | -5 to 105 °C, 23 to 221 °F  |
|             | Resolution              | 0.1   |
| Temperature | Relative<br>Accuracy    | ±0.1  |
|             | Offset<br>Calibration   | 1 point   |
| Datalogging | Number of Datalogs      | 5000 with time and date stamp   |
|             | Log Function            | Manual, ready (include AUTO-READ $^{\text{\tiny{TM}}}$ ), timed   |
| Inputs      | pH Electrode            | BNC, reference pin  |
| mpato       | ATC Probe               | 8-pin mini-DIN  |
| Output      |                         | RS232, USB  |
|             | AC Adapter              | Optional – universal, 100-240 VAC   |
| Power       | Battery<br>Power        | Included – 4 AAs  |
|             | Battery Life            | 800 hrs   |
|             |                         |   |

## **Meter Ordering Information**

| Cat. No.  | Description  |
|-----------|--|
| STARA3240 | Orion Star A324 pH/ISE Portable Meter - Four AA batteries (installed)  |
| STARA3245 | Orion Star A324 pH/ISE Portable Meter Kit  - 8107UWMMD Orion ROSS Triode 3-in-1 pH/ATC Probe - Orion pH 4, 7 and 10 buffers, 10 pouches of each value - 810001 Orion ROSS electrode storage solution - 911110 Rinse solution, 10 pouches - Protective armor with electrode holder - Hard carrying case - Four AA batteries (installed) |

## Accessories

| Cat. No.   | Description  |
|------------|--|
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |
| 1210005    | Soft carrying case for Orion Star, Star Plus and Star A-series portable meters                               |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star<br>A-series armor  |
| 8107UWMMD  | Orion ROSS Triode 3-in-1 pH/ATC probe, gel-filled, 3m cable  |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |
| 910410     | Orion pH 4.01 buffer pouches, 10/pk  |
| 910710     | Orion pH 7.00 buffer pouches, 10/pk  |
| 911010     | Orion pH 10.01 buffer pouches, 10/pk   |
| 911110     | Orion pH rinse solution pouches, 10/pk   |
| 9512HPBNWP | Orion high-performance ammonia ISE   |
| 9512BNWP   | Orion standard ammonia ISE   |
| 951007     | Ammonia standard, 1000 ppm, 475 mL   |
| 951210     | Low-level ammonia ISA  |
| 951211     | Ammonia ISA, 475 mL  |
| 951209     | Fill solution for Orion high-performance ammonia ISE, 60 mL  |
| 951202     | Fill Solution for standard ammonia ISE, 60 mL  |
| 951213     | Ammonia ISE storage solution   |
| 9609BNWP   | Orion fluoride combination ISE   |
| 040906     | Fluoride standard, 1 ppm with TISAB II, 475 mL   |
| 040907     | Fluoride standard, 2 ppm with TISAB II, 475 mL   |
| 040908     | Fluoride standard, 10 ppm with TISAB II, 475 mL  |
| 940909     | TISAB II, 1 gallon   |
| 8611BNWP   | Orion ROSS sodium combination ISE  |

рН 📕

ISE

ORP

Cond

RDO 🔣

DO 📗





pH/Conductivity Portable Meter

The Thermo Scientific Orion Star A325 pH/Conductivity Portable Meter is the top choice for advanced pH and conductivity measurements on location. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows individual pH/mV/RmV or conductivity/TDS/salinity/resistivity results with temperature, or all results at once. Additional helpful information include electrode status, time, date, sample ID, user ID and calibration points. Advanced features such as stability and averaging options provide additional options allowing for laboratory accuracy everywhere you need it.

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Fix pH calibration errors without a complete recalibration with calibration editing
- Selectable cell constant allows for use with 2- or 4-cell conductivity cells
- For advanced accuracy in conductivity readings, selectable reading reference temperatures of 5, 10, 15, 20 or 25 °C with linear, non-linear, nLFu or USP/EP curve options
- Choice of practical salinity or natural sea water curves for salinity readings
- Linear or ISO/EN 27888 curves for TDS
- Non-volatile memory holds up to 5000 data points with time and date stamp
- User and sample ID fields for additional data accuracy
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty





| Product Sp   | ecifications                                     | <b>i</b>                                     |
|--------------|--|--|
|              |  | CE   TUV   FCC                               |
|              | Range  | -2.000 to 20.000                             |
| рН           | Resolution                                       | 0.1, 0.01, 0.001                             |
|              | Relative<br>Accuracy                             | ±0.002                                       |
|              | Calibration<br>Points                            | Up to 5                                      |
|              | Calibration<br>Editing                           | Yes  |
|              | Range  | ±2000.0 mV                                   |
| mV/Rel mV/   | Resolution                                       | 0.1  |
| ORP          | Relative   | ±0.2 mV or ±0.05 % of reading                |
|              | Accuracy   | whichever is greater                         |
|              | E <sub>H</sub> ORP Mode<br>Range                 | Yes<br>0.001 µS to 3000 mS                   |
|              |  | 0.001 µS minimum;                            |
|              | Resolution                                       | 4 significant figures minimum                |
|              | Relative<br>Accuracy                             | 0.5 % reading ±1 digit                       |
| Conductivity | Reference<br>Temperature                         | 5, 10, 15, 20, 25 °C (default)               |
| Conductivity | Temperature<br>Compensation                      | Linear (0 to 10.0 %/°C), nLn, nLFu, USP/EP   |
|              | Compatible<br>Cell Constants                     | 0.001 to 199.9                               |
|              | Calibration<br>Points                            | Up to 5 points                               |
|              | Calibration<br>Editing                           | Yes  |
|              | Range  | 2 ohm to 100 meg-ohm                         |
| Resistivity  | Resolution                                       | 2 ohms-cm                                    |
|              | Relative<br>Accuracy                             | 0.5 % reading ±1 digit                       |
|              | Туре   | Practical salinity or natural sea water      |
|              | Range  | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water |
| Salinity     | Resolution                                       | 0.01   |
|              | Relative<br>Accuracy                             | ±0.1   |
|              | Range  | 0 to 200 ppt                                 |
|              | Resolution                                       | 4 significant digits                         |
| TDS          | Relative<br>Accuracy                             | ±0.5 % reading ±1 digit                      |
|              | TDS Factor<br>Range                              | Linear 0.01 to 10.00, default 0.49           |
|              | Range  | -5 to 105 °C, 23 to 221 °F                   |
|              | Resolution                                       | 0.1  |
| Temperature  | Relative<br>Accuracy                             | ±0.1   |
|              | Offset<br>Calibration                            | 1 point                                      |
| Datalogging  | Number of Datalogs                               | 5000 with time and date stamp                |
|              | Log Function                                     | Manual, ready (include AUTO-READ™), timed    |
|              | pH Electrode                                     | BNC, reference pin                           |
| Inputs       | Conductivity with Temperature Probe or ATC Probe | 8-pin mini-DIN                               |
| Output       |  | RS232, USB                                   |
|              | AC Adapter                                       | Optional – universal, 100-240 VAC            |
| Power        | Battery Power                                    | Included – 4 AAs                             |
|              | Battery Life                                     | 800 hrs                                      |
|              |  |  |

## **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA3250 | Orion Star A325 pH/Conductivity Portable Multiparameter Meter - Four AA batteries (installed)   |
| STARA3255 | Orion Star A325 pH/Conductivity Portable Multiparameter Meter Kit  - 8107UWMMD Orion ROSS Triode 3-in-1 pH/ATC probe, 3m cable - 013010MD Orion DuraProbe conductivity cell, 3m cable - 910410 Orion pH 4.01 buffer, 10 pouches - 910710 Orion pH 7.00 buffer, 10 pouches - 911010 Orion pH 10.01 buffer, 10 pouches - 911010 Orion pH 10.01 buffer, 10 pouches - 810001 Orion ROSS electrode storage solution, 475 mL - 01100710 Orion 1413 µS conductivity standard, 10 pouches - 911110 rinse solution, 10 pouches - Protective armor with electrode holder - Hard carrying case - Four AA batteries (installed) |

#### Accessories

| Cat. No.   | Description  |  |  |
|------------|--|--|--|
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |  |  |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |  |  |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |  |  |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star<br>A-series armor  |  |  |
| 8107UWMMD  | Orion ROSS Triode 3-in-1 pH/ATC probe, gel-filled,<br>3m cable   |  |  |
| 013010MD   | Orion DuraProbe conductivity cell with 3m cable, 4-cell, K=0.475   |  |  |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |  |  |
| 910410     | Orion pH 4.01 buffer pouches, 10/pk  |  |  |
| 910710     | Orion pH 7.00 buffer pouches, 10/pk  |  |  |
| 911010     | Orion pH 10.01 buffer pouches, 10/pk   |  |  |
| 911110     | Orion pH rinse solution pouches, 10/pk   |  |  |
| 011008     | Orion conductivity standard, 100 µS/cm, 5 x 60 mL bottles  |  |  |
| 01100710   | Orion conductivity standard pouches, 1413 µS/cm, 10/pk   |  |  |
| 01100610   | Orion conductivity standard pouches, 12.9 mS/cm, 10/pk   |  |  |
| 1010001    | Orion conductivity calibration resistor kit  |  |  |



рН

ISE

ORP

Cond

RDO





pH/RDO/Dissolved Oxygen Portable Meter

The Thermo Scientific Orion Star A326 pH/RDO/Dissolved Oxygen Portable Meter is the top choice for advanced pH and dissolved oxygen measurements on location. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows individual pH/mV/RmV or dissolved oxygen in percent saturation or concentration results with temperature, or all results at once. Additional helpful information include electrode status, time, date, sample ID, user ID and calibration points. Advanced features such as stability and averaging options provide additional options allowing for laboratory accuracy everywhere you need it.



- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Fix pH calibration errors without a complete recalibration with calibration editing
- No need to worry about DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- Easily get accurate DO results with automatic temperature, manual salinity and automatic/manual barometric pressure correction
- Calibrate DO using water-saturated air, air-saturated water, a Winkler titration or zero point calibration
- Non-volatile memory holds up to 5000 data points with time and date stamp
- User and sample ID fields for additional data accuracy
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use on AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty



| Product Sp          | ecifications                                   | •   |
|---------------------|--|---|
|                     |  | CE   TUV   FCC  |
| рН                  | Range  | -2.000 to 20.000  |
|                     | Resolution                                     | 0.1, 0.01, 0.001  |
|                     | Relative<br>Accuracy                           | ±0.002  |
|                     | Calibration<br>Points                          | Up to 5   |
|                     | Calibration<br>Editing                         | Yes   |
|                     | Range  | ±2000.0 mV  |
|                     | Resolution                                     | 0.1   |
| mV/Rel mV/<br>ORP   | Relative                                       | ±0.2 mV or ±0.05 % of reading   |
|                     | Accuracy                                       | whichever is greater  |
|                     | E <sub>H</sub> ORP Mode                        | Yes   |
|                     | _  | - Polarographic   |
|                     | Range  | 0 to 90 mg/L  |
|                     | Resolution                                     | 0.01, 0.1   |
|                     | Relative<br>Accuracy                           | ±0.2  |
|                     | _  | Polarographic   |
|                     | Range  | 0 to 600  |
|                     | Resolution                                     | 0.1, 1  |
|                     | Relative<br>Accuracy                           | ±2 %  |
|                     | Concentration                                  | - RDO   |
|                     | Range  | 0 to 50 mg/L  |
|                     | Resolution                                     | 0.01, 0.1   |
| Dissolved<br>Oxygen | Relative<br>Accuracy                           | ±0.1 mg/L up to 8 mg/L;<br>±0.2 mg/L from 8 to 20 mg/L;<br>10 % of reading from 20 to 50 mg/L |
| охуусп              | % Saturation -                                 |   |
|                     | Range  | 0 to 500  |
|                     | Resolution                                     | 0.1, 1  |
|                     | Relative<br>Accuracy                           | ±2 %  |
|                     | Manual<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg  |
|                     | Manual<br>Salinity Factor<br>Correction        | 0 to 45 ppt   |
|                     | Calibration<br>Features                        | Water-saturated air, air saturated water,<br>manual (Winkler) and zero point                  |
|                     | Probe<br>Characteristics                       | Polarographic or RDO  |
|                     | Range  | 0 to 50 °C, 32 to 122 °F  |
|                     | Resolution                                     | 0.1   |
| Temperature         | Relative<br>Accuracy                           | ±0.1  |
|                     | Offset<br>Calibration                          | 1 point   |
| Datalogging         | Number of Datalogs                             | 5000 with time and date stamp   |
| 353                 | Log Function                                   | Manual, ready (includes AUTO-READ $^{\mathrm{m}}$ ), timed                                    |
|                     | pH Electrode                                   | BNC, reference pin  |
| Inputs              | DO with<br>Temp. Probe                         | 9-pin mini-DIN  |
| Output              |  | RS232, USB  |
|                     | AC Adapter                                     | Optional – universal, 100-240 VAC   |
| Power               | <b>Battery Power</b>                           | Optional – 4 AAs  |
|                     | Battery Life                                   | 800 hrs   |
|                     |  |   |

## **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA3260 | Orion Star A326 pH/RDO/Dissolved Oxygen Portable<br>Multiparameter Meter<br>- Four AA batteries (installed)   |
| STARA3265 | Orion Star A326 pH/RD0/Dissolved Oxygen Portable Multiparameter Meter Kit  - 8107UWMMD Orion ROSS Triode 3-in-1 pH/ATC probe, 3m cable  - 087010MD Orion RDO optical/luminescence-based probe with 3m cable, stainless steel guard, calibration sleeve and optical/luminescence-based cap  - 910410 Orion pH 4.01 buffer, 10 pouches  - 910710 Orion pH 7.00 buffer, 10 pouches  - 911010 Orion pH 10.01 buffer, 10 pouches  - 810001 Orion ROSS electrode storage solution, 475 mL  - 01100710 Orion 1413 µS conductivity standard, 10 pouches  - 911110 rinse solution, 10 pouches  - Protective armor with electrode holder  - Hard carrying case  - Four AA batteries (installed) |

## Accessories

|            | ~  |
|------------|--|
| Cat. No.   | Description  |
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star<br>A-series armor  |
| 8107UWMMD  | Orion ROSS Triode 3-in-1 pH/ATC probe, gel-filled, 3m cable  |
| 087020MD   | Orion RDO optical DO probe, 6m cable with stainless steel guard and calibration sleeve                       |
| 087030MD   | Orion RDO optical DO probe, 10m cable with stainless steel guard and calibration sleeve                      |
| 087050MD   | Orion RDO optical DO probe, 15m cable with stainless steel guard and calibration sleeve                      |
| 087100MD   | Orion RDO optical DO probe, 30m cable with stainless steel guard and calibration sleeve                      |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |
| 910410     | Orion pH 4.01 buffer pouches, 10/pk  |
| 910710     | Orion pH 7.00 buffer pouches, 10/pk  |
| 911010     | Orion pH 10.01 buffer pouches, 10/pk   |
| 911110     | Orion pH rinse solution pouches, 10/pk   |
| 087010MD   | Orion RDO optical DO probe, 3m cable with stainless steel guard and calibration sleeve                       |
| 083005MD   | Orion polarographic DO probe, 1.5m cable   |









Cond

RDO 📉

DO 🗾





# Thermo Scientific Orion Star A321 & Star A221

pH Portable Meters

## Orion Star A321

The Thermo Scientific Orion Star A321 pH Portable Meter is the top choice for advanced pH measurements on location. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows pH/mV/RmV and temperature along with electrode status, time, date, sample ID, user ID and calibration points. Advanced features such as stability and averaging options provide additional options allowing for laboratory accuracy everywhere you need it.

## **Features and Benefits**

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Fix errors without a complete recalibration with calibration editing
- Non-volatile memory holds up to 5000 data points with time and date stamp
- User and sample ID fields for additional data accuracy
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty



#### Orion Star A221

Get the readings that you need where you need them with the Thermo Scientific Orion Star A221 pH Portable Meter. The large, backlit graphic LCD display shows pH/mV/RmV and temperature along with electrode status, time, date and calibration points.

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers
- Non-volatile memory holds up to 1000 data points with time and date stamp
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty





## Product Specifications (Orion Star A321 & Star A221)

| Comodione               | (Orion Star A321 & Star A221)  |
|-------------------------|--|
|                         | CE   TUV   FCC   |
| Range                   | -2.000 to 20.000   |
| Resolution              | 0.1, 0.01, 0.001   |
| Relative<br>Accuracy    | ±0.002   |
| Calibration<br>Points   | Up to 5  |
| Calibration<br>Editing  | Orion Star A321: Yes<br>Orion Star A221: No  |
| Range                   | ±2000.0 mV   |
| Resolution              | 0.1  |
| Relative<br>Accuracy    | ±0.2 mV or ±0.05 % of reading whichever is greater   |
| E <sub>H</sub> ORP Mode | Yes  |
| Range                   | -5 to 105 °C, 23 to 221 °F   |
| Resolution              | 0.1  |
| Relative<br>Accuracy    | ±0.1   |
| Offset<br>Calibration   | 1 point  |
| Number<br>of Datalogs   | Orion Star A321:<br>5000 with time and date stamp<br>Orion Star A221:<br>1000 with time and date stamp   |
| Log Function            | Orion Star A321:<br>Manual, ready (includes AUTO-READ™), timed<br>Orion Star A221:<br>Manual, Automatic with AUTO-READ™  |
| pH Electrode            | BNC, reference pin   |
| ATC Probe               | 8-pin mini-DIN   |
|                         | RS232, USB   |
| AC Adapter              | Optional – universal, 100-240 VAC  |
| Battery<br>Power        | Included – 4 AAs   |
| <b>Battery Life</b>     | 800 hrs  |
|                         | Range Resolution Relative Accuracy Calibration Points Calibration Editing Range Resolution Relative Accuracy E <sub>H</sub> ORP Mode Range Resolution Relative Accuracy Offset Calibration Number of Datalogs  Log Function pH Electrode ATC Probe  AC Adapter Battery Power |

## **Meter Ordering Information**

| Cat. No.  | Description  |
|-----------|--|
| STARA3210 | Orion Star A321 pH Portable Meter - Four AA batteries (installed)  |
| STARA3215 | Orion Star A321 pH Portable Meter Kit  - 8107UWMMD Orion ROSS Triode 3-in-1 pH/ATC probe - 910410 Orion pH 4.01 buffer, 10 pouches - 910710 Orion pH 7.00 buffer, 10 pouches - 911010 Orion pH 10.01 buffer, 10 pouches - 810001 Orion ROSS electrode storage solution, 475 mL - 911110 rinse solution, 10 pouches - Protective armor with electrode holder - Hard carrying case - Four AA batteries (installed) |
| STARA2210 | Orion Star A221 pH Portable Meter - Four AA batteries (installed)  |
| STARA2215 | Orion Star A221 pH Portable Meter Kit  - 8107UWMMD Orion ROSS Triode 3-in-1 pH/ATC probe - 910410 Orion pH 4.01 buffer, 10 pouches - 910710 Orion pH 7.00 buffer, 10 pouches - 911010 Orion pH 10.01 buffer, 10 pouches - 810001 Orion ROSS electrode storage solution, 475 mL - 911110 rinse solution, 10 pouches - Protective armor with electrode holder - Hard carrying case - Four AA batteries (installed) |

## Accessories

| Accessories  |  |  |
|--|--|--|
| Description  |  |  |
| Hard carrying case for Orion Star A-series portable meters   |  |  |
| Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |  |  |
| pH electrode holder for Orion Star A-series armor  |  |  |
| Conductivity and DO probe holder for Orion Star<br>A-series armor  |  |  |
| Orion ROSS Triode 3-in-1 pH/ATC probe, gel-filled, 3m cable  |  |  |
| Optional universal power adapter, 100-240 V, 50/60 Hz  |  |  |
| Orion pH 4.01 buffer pouches, 10/pk  |  |  |
| Orion pH 7.00 buffer pouches, 10/pk  |  |  |
| Orion pH 10.01 buffer pouches, 10/pk   |  |  |
| Orion pH rinse solution pouches, 10/pk   |  |  |
|  |  |  |





ISE \_

\_\_\_

ORP

Cond

RDO 📉





pH Portable Meter

Thermo Scientific Orion Star A121 pH Portable Meters combine simplicity with accuracy. A large LCD displays pH or mV readings along with temperature. Icons provide quick updates on battery life, electrode status and calibration information. Simple button layout and onscreen messages help with calibration and setup menu choices.



- Large, informative screen clearly shows the main information, including pH or mV reading, temperature in °C or °F, meter mode, calibration information and battery life
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Up to 3 point pH calibration
- Automatic recognition of USA/NIST and DIN buffers
- Easily recall calibration slope data for procedure checks
- Non-volatile memory holds up to 50 data points
- Four AA batteries (included) provide over 2000 hours of operation or purchase the universal power adapter (sold separately) to run meter on AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty



|             |                                      | CE   TUV   FCC                                     |
|-------------|--------------------------------------|--|
|             | Range                                | -2.00 to 16.00                                     |
| рН          | Resolution                           | 0.1, 0.01  |
|             | Relative<br>Accuracy                 | ±0.01  |
|             | pH<br>Calibration<br>Points          | Up to 3  |
|             | Range - mV                           | ±1600.0 mV   |
|             | Range - RmV                          | ±1999.9 mV   |
| mV/Rel mV/  | Resolution                           | 0.1  |
| ORP         | Relative<br>Accuracy                 | ±0.2 mV or ±0.05 % of reading whichever is greater |
|             | E <sub>H</sub> ORP Mode              | No   |
|             | Range                                | -5 to 105 °C, 23 to 221 °F                         |
|             | Resolution                           | 0.1  |
|             | Relative<br>Accuracy                 | ±0.1   |
| Temperature | Temperature<br>Compensation          | Automatic or manual                                |
|             | Temperature<br>Offset<br>Calibration | Yes  |
| Datalogging | Number of Datalogs                   | 50   |
|             | Log Function                         | Manual, Automatic with AUTO-READ™                  |
| Inputs      | pH Electrode                         | BNC  |
| IIIputo     | ATC Probe                            | 8-pin mini-DIN                                     |
|             | AC Adapter                           | Optional – universal, 100-240 VAC                  |
| Power       | Battery<br>Power                     | Included – 4 AAs                                   |
|             | Battery Life                         | 2000 hrs   |

## **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA1210 | Orion Star A121 pH Portable Meter<br>- Four AA batteries (installed)  |
| STARA1215 | Orion Star A121 pH Portable Meter Kit - 9107BNMD Orion Triode 3-in-1 pH/ATC probe - Solution kit (60 mL each of pH 4, 7 and 10 buffers; storage solution and cleaning solution) - 911110 rinse solution, 10 pouches - Protective armor with electrode holder - Hard carrying case - Four AA batteries (installed) |

## Accessories

| Cat. No.   | Description  |
|------------|--|
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes     |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star<br>A-series armor  |
| 9157BNMD   | Orion Triode 3-in-1 pH/ATC probe, refillable, epoxy body   |
| 9107BNMD   | Orion Triode 3-in-1 pH/ATC probe, gel-filled, epoxy body   |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |
| 910410     | Orion pH 4.01 buffer pouches, 10/pk  |
| 910710     | Orion pH 7.00 buffer pouches, 10/pk  |
| 911010     | Orion pH 10.01 buffer pouches, 10/pk   |
| 911110     | Orion pH rinse solution pouches, 10/pk   |
| 916099     | Orion pH buffer assortment pack - 60 mL each of pH 4, 7 and 10 buffers; storage solution and cleaning solution D |

рН 📕

ISE \_

ORP

Cond

RDO 📗

DO





# Thermo Scientific Orion Star A322 & Star A222

Conductivity Portable Meters

## Orion Star A322

The Thermo Scientific Orion Star A322 Conductivity Portable Meter is the top choice for advanced measurements on location. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows conductivity/ TDS/salinity resistivity and temperature along with electrode status, time, date, sample ID, user ID and calibration points. Advanced features such as stability and averaging options provide additional options allowing for laboratory accuracy everywhere you need it.

## **Features and Benefits**

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Selectable cell constant allows for use with 2- or 4-cell conductivity cells
- For advanced accuracy, selectable reading reference temperatures of 5, 10, 15, 20 or 25 °C with linear, non-linear, nLFu or USP/EP curve options
- Choice of practical salinity or natural sea water curves for salinity readings
- Linear or ISO/EN 27888 curves for TDS
- Non-volatile memory holds up to 5000 data points with time and date stamp
- User and sample ID fields for additional data accuracy
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty

## Orion Star A222

Get the readings that you need where you need them with the Thermo Scientific Orion Star A222 Conductivity Portable Meter. The large, backlit graphic LCD display shows conductivity/TDS/salinity/resistivity and temperature along with electrode status, time, date and calibration points.

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Selectable cell constant allows for use with 2- or 4-cell conductivity cells
- Selectable reading reference temperatures of 15, 20 or 25 °C with linear or non-linear curve options
- Includes the standard practical salinity curve for salinity readings
- Linear or ISO/EN 27888 curves for TDS
- Non-volatile memory holds up to 1000 data points with time and date stamp
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty







## Product Specifications (Orion Star A322 & Star A222)

| Product Sp   | ecifications                                 | (Orion Star A322 & Star A222)   |
|--------------|--|---|
|              |  | CE   TUV   FCC  |
|              | Range  | 0.001 µS to 3000 mS   |
|              | Resolution                                   | 0.001 µS minimum;<br>4 significant figures minimum  |
|              | Relative<br>Accuracy                         | 0.5 % reading ±1 digit  |
|              | Reference<br>Temperature                     | Orion Star A322: 5, 10, 15, 20, 25 °C<br>Orion Star A222: 15, 20, 25 °C   |
| Conductivity | Temperature<br>Compensation                  | Orion Star A322:<br>Linear (0 to 10.0 %/°C), nLn, nLFu, USP/EP<br>Orion Star A222:<br>Linear (0 to 10.0 %/°C), nLn      |
|              | Compatible<br>Cell Constants                 | 0.001 to 199.9  |
|              | Calibration<br>Points                        | Up to 5 points  |
|              | Calibration<br>Editing                       | Orion Star A322: Yes<br>Orion Star A222: –  |
|              | Range  | 2 ohm to 100 meg-ohm  |
| Resistivity  | Resolution                                   | 2 ohms-cm   |
|              | Relative<br>Accuracy                         | 0.5 % reading ±1 digit  |
|              | Туре   | Orion Star A322:<br>Practical salinity or natural sea water<br>Orion Star A222:<br>Practical salinity                   |
| Salinity     | Range  | 0.01 to 80.0 psu<br>0.01 to 42 ppt sea water  |
|              | Resolution                                   | 0.01  |
|              | Relative<br>Accuracy                         | ±0.1  |
|              | Range  | 0 to 200 ppt  |
|              | Resolution                                   | 4 significant digits  |
| TDS          | Relative<br>Accuracy                         | ±0.5 % reading ±1 digit   |
|              | TDS Factor<br>Range                          | Linear, 0.01 to 10.00   |
|              | Range  | -5 to 105 °C, 23 to 221 °F  |
|              | Resolution                                   | 0.1   |
| Temperature  | Relative<br>Accuracy                         | ±0.1  |
|              | Offset<br>Calibration                        | 1 point   |
| Datalogging  | Number<br>of Datalogs                        | Orion Star A322:<br>5000 with time and date stamp<br>Orion Star A222:<br>1000 with time and date stamp                  |
|              | Log Function                                 | Orion Star A322:<br>Manual, ready (includes AUTO-READ™), timed<br>Orion Star A222:<br>Manual, Automatic with AUTO-READ™ |
| Inputs       | Conductivity<br>with<br>Temperature<br>Probe | 8-pin mini-DIN  |
| Output       |  | RS232, USB  |
|              | AC Adapter                                   | Optional – universal, 100-240 VAC   |
| Power        | <b>Battery Power</b>                         | Included – 4 AAs  |
|              | Battery Life                                 | 800 hrs   |

## **Meter Ordering Information**

| Cat. No.  | Description   |
|-----------|---|
| STARA3220 | Orion Star A322 Conductivity Portable Meter - Four AA batteries (installed)   |
| STARA3225 | Orion Star A322 Conductivity Portable Meter Kit - 013010MD Orion DuraProbe conductivity cell with 3m cable, K=0.475 - 01100710 conductivity standard pouches, 1413 μS, pack of 10 - 911110 rinse solution, pack of 10 - Protective armor with electrode holder - Hard carrying case - Four AA batteries (installed)       |
| STARA2220 | Orion Star A222 Conductivity Portable Meter - Four AA batteries (installed)   |
| STARA2225 | Orion Star A222 Conductivity Portable Meter Kit  - 013010MD Orion DuraProbe conductivity cell with 3m cable, K=0.475  - 01100710 conductivity standard pouches, 1413 μS, pack of 10  - 911110 rinse solution, pack of 10  - Protective armor with electrode holder  - Hard carrying case  - Four AA batteries (installed) |

## **Accessories**

| Cat. No.   | Description  |
|------------|--|
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |
| 1210005    | Soft carrying case for Orion Star, Star Plus and Star A-series portable meters                               |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star<br>A-series armor  |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |
| 013005MD   | Orion DuraProbe conductivity cell with 1.5m cable, 4-cell, K=0.475   |
| 013010MD   | Orion DuraProbe conductivity cell with 3m cable, 4-cell, K=0.475   |
| 013020MD   | Orion DuraProbe conductivity cell with 6m cable, 4-cell, K=0.475   |
| 013025MD   | Orion DuraProbe conductivity cell with 10m cable, 4-cell, K=0.475  |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |
| 011008     | Orion conductivity standard, 100 $\mu$ S/cm, 5 x 60 mL bottles   |
| 01100710   | Orion conductivity standard pouches, 1413 µS/cm, 10/pk   |
| 01100610   | Orion conductivity standard pouches, 12.9 mS/cm, 10/pk   |
| 911110     | Orion rinse solution pouches, 10/pk  |
| 1010001    | Orion conductivity calibration resistor kit  |













RDO 📗







Conductivity Portable Meter

Thermo Scientific Orion Star A122 Conductivity Portable Meters combine simplicity with accuracy. A large LCD displays conductivity or TDS readings along with temperature. Icons provide quick updates on battery life and completed calibration status. Simple button layout and onscreen messages help with calibration and setup menu choices.

- Large, informative screen clearly shows the main information, including conductivity or TDS reading, temperature in °C or °F, meter mode and battery life
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Selectable cell constant
- Manual calibration or automatic calibration when using a Thermo Scientific Orion conductivity standard (sold separately)
- Selectable reading reference temperatures of 20 °C or 25 °C for an accurate result
- Non-volatile memory holds up to 50 data points
- Four AA batteries (included) provide over 2000 hours of operation or purchase the universal power adapter (sold separately) to run meter on AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty





| i roduct op  | comeations                   |   |
|--------------|------------------------------|---|
|              |                              | CE   TUV   FCC                                    |
|              | Range                        | 0.1 µS/cm to 200 mS/cm                            |
|              | Resolution                   | 0.01 µS minimum;<br>3 significant figures minimum |
|              | Relative<br>Accuracy         | 0.5 % reading ±1 digit                            |
| Conductivity | Reference<br>Temperature     | 20, 25 °C (default)                               |
|              | Temperature<br>Compensation  | Linear (0.0 to 10.0 %/°C)                         |
|              | Compatible<br>Cell Constants | 0.001 to 10                                       |
|              | Calibration<br>Points        | 1   |
|              | Range                        | 1 to 19,999 mg/L                                  |
|              | Resolution                   | 4 significant digits                              |
| TDS          | Relative<br>Accuracy         | 0.5 % reading ±1 digit                            |
|              | TDS Factor<br>Range          | Linear, 0.01 to 10.00, default 0.49               |
|              | Range                        | -5 to 105 °C, 23 to 221 °F                        |
|              | Resolution                   | 0.1   |
| Temperature  | Relative<br>Accuracy         | ±0.1  |
|              | Offset<br>Calibration        | Yes   |
| Datalogging  | Number of Datalogs           | 50  |
|              | Log Function                 | Manual, Automatic with AUTO-READ™                 |
| Inputs       | Conductivity                 | 8-pin mini-DIN                                    |
|              | AC Adapter                   | Optional – universal, 100-240 VAC                 |
| Power        | <b>Battery Power</b>         | Included – 4 AAs                                  |
|              | Battery Life                 | 2000 hrs  |
|              |                              |   |

## **Meter Ordering Information**

| Cat. No.  | Description  |
|-----------|--|
| STARA1220 | Orion Star A122 Conductivity Portable Meter - Four AA batteries (installed)  |
| STARA1225 | Orion Star A122 Conductivity Portable Meter Kit - 011050MD Orion conductivity cell, K=1.0 - 011007 conductivity standard, 1413 µS, 5 x 60 mL - 911110 rinse solution, 10 pouches - Protective armor with electrode holder - Hard carrying case - Four AA batteries (installed) |

## Accessories

| Cat. No.   | Description  |  |  |  |
|------------|--|--|--|--|
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |  |  |  |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |  |  |  |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star<br>A-series armor  |  |  |  |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |  |  |  |
| 011050MD   | Orion 2-electrode conductivity cell, K=1.0   |  |  |  |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |  |  |  |
| 011008     | Orion conductivity standard, 100 µS/cm, 5 x 60 mL bottles  |  |  |  |
| 01100710   | Orion conductivity standard pouches, 1413 µS/cm, 10/pk   |  |  |  |
| 01100610   | Orion conductivity standard pouches, 12.9 mS/cm, 10/pk   |  |  |  |
| 911110     | Rinse solution, 10/pk  |  |  |  |
| 1010001    | Orion conductivity calibration resistor kit  |  |  |  |

рН 📗

ISF

ORP |

Cond

RDO

DO 📗





## Thermo Scientific Orion Star A323 & Star A223

RDO/Dissolved Oxygen Portable Meters

## Orion Star A323

The Thermo Scientific Orion Star A323 RDO/Dissolved Oxygen Portable Meter is the top choice for ensuring proper oxygen levels for maximum process control on location. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows dissolved oxygen in percent saturation or concentration with temperature and electrode status, time, date, sample ID, user ID and calibration points. Advanced features such as stability and averaging options provide additional options allowing for laboratory accuracy everywhere you need it.

#### **Features and Benefits**

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- No need to worry about DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- Accurate results are easy with automatic temperature, manual salinity and automatic/manual barometric pressure compensation
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen, ready indicator alerts when readings are stable and timed reading gathers data in specific time intervals
- Calibrate using water-saturated air, air-saturated water, a Winkler titration or zero point calibration
- Non-volatile memory holds up to 5000 data points with time and date stamp
- User and sample ID fields for additional data accuracy
- Stores up to 10 methods allow for customized procedures to differentiate between tests and/or users
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty



#### Orion Star A223

Get the readings that you need where you need them with the Thermo Scientific Orion Star A223 RDO/Dissolved Oxygen Portable Meter. The large, backlit graphic LCD display shows dissolved oxygen in percent saturation or concentration with temperature and electrode status, time, date and calibration points.

- Anyone can use this meter with its plain language screen prompts, soft keys that update for easy selection and a multilanguage interface that supports English, Spanish, French, Italian, German and Chinese
- No need to worry about DO probe compatibility the meter accepts and automatically recognizes Orion polarographic and RDO optical/luminescence-based probes
- Accurate results are easy with automatic temperature, manual salinity and automatic/manual barometric pressure compensation
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Calibrate using water-saturated air, air-saturated water, a Winkler titration or zero point calibration
- Non-volatile memory holds up to 1000 data points with time and date stamp
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary data analysis software
- Four AA batteries (included) provide over 800 hours of operation or purchase the universal power adapter (sold separately) to use AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty





| Product Specifications (Orion Star A323 & Star A223) |  |   |  |  |
|--|--|---|--|--|
|  |  | CE   TUV   FCC  |  |  |
|  | Concentration                                  | - Polarographic   |  |  |
|  | Range  | 0 to 90 mg/L  |  |  |
|  | Resolution                                     | 0.01, 0.1   |  |  |
|  | Relative<br>Accuracy                           | ±0.2  |  |  |
|  | % Saturation -                                 | Polarographic   |  |  |
|  | Range  | 0 to 600  |  |  |
|  | Resolution                                     | 0.1, 1  |  |  |
|  | Relative<br>Accuracy                           | ±2 %  |  |  |
|  | Concentration                                  | - RDO   |  |  |
|  | Range  | 0 to 50 mg/L  |  |  |
|  | Resolution                                     | 0.01, 0.1   |  |  |
| Dissolved<br>Oxygen                                  | Relative<br>Accuracy                           | ±0.1 mg/L up to 8 mg/L;<br>±0.2 mg/L from 8 to 20 mg/L;<br>10 % of reading from 20 to 50 mg/L                             |  |  |
|  | % Saturation -                                 | · RDO   |  |  |
|  | Range  | 0 to 500  |  |  |
|  | Resolution                                     | 0.1, 1  |  |  |
|  | Relative<br>Accuracy                           | ±2 %  |  |  |
|  | Manual<br>Barometric<br>Pressure<br>Correction | 450.0 to 850.0 mm Hg  |  |  |
|  | Manual<br>Salinity Factor<br>Correction        | 0 to 45 ppt   |  |  |
|  | Calibration<br>Features                        | Water saturated air, air saturated water, manual (Winkler) and zero point   |  |  |
|  | Probe<br>Characteristics                       | Polarographic or RDO  |  |  |
|  | Range  | 0 to 50 °C, 32 to 122 °F  |  |  |
|  | Resolution                                     | 0.1   |  |  |
| Temperature  | Relative<br>Accuracy                           | ±0.1  |  |  |
|  | Offset<br>Calibration                          | 1 point   |  |  |
| Datalogging  | Number of<br>Datalogs                          | Orion Star A323:<br>5000 with time and date stamp<br>Orion Star A223:<br>1000 with time and date stamp                    |  |  |
| Dutilogging  | Log Function                                   | Orion Star A323:<br>Manual, automatic (with Auto-Read™), timed<br>Orion Star A223:<br>Manual, automatic (with Auto-Read™) |  |  |
| Inputs   | DO with<br>Temperature<br>Probe                | 9-pin mini-DIN  |  |  |
| Output   |  | RS232, USB  |  |  |
|  | <b>Power Cord</b>                              | Optional – universal, 100-240 VAC   |  |  |
| Power  | Battery<br>Power                               | Included – 4 AAs  |  |  |
|  | <b>Battery Life</b>                            | 800 hrs   |  |  |

## **Meter Ordering Information**

| 3         |  |  |  |  |
|-----------|--|--|--|--|
| Cat. No.  | Description  |  |  |  |
| STARA3230 | Orion Star A323 Dissolved Oxygen Portable Meter - Four AA batteries (installed)  |  |  |  |
| STARA3235 | Orion Star A323 Dissolved Oxygen Portable Meter Kit  - 087010MD Orion RDO optical/luminescence-based probe with 3m cable, stainless steel guard, calibration sleeve and optical/luminescence-based cap  - Protective armor with electrode holder  - Hard carrying case  - Four AA batteries (installed)                |  |  |  |
| STARA2230 | Orion Star A223 RDO/Dissolved Oxygen Portable Meter - Four AA batteries (installed)  |  |  |  |
| STARA2235 | Orion Star A223 RDO/Dissolved Oxygen Portable Meter Kit  - 087010MD Orion RDO optical/luminescence-<br>based probe with 3m cable, stainless<br>steel guard, calibration sleeve and optical/<br>luminescence-based cap  - Protective armor with electrode holder  - Hard carrying case  - Four AA batteries (installed) |  |  |  |

## Accessories

| Accessories   |  |  |  |  |
|---|--|--|--|--|
| Cat. No.  | Description  |  |  |  |
| STARA-CS  | Hard carrying case for Orion Star A-series portable meters   |  |  |  |
| STARA-AR  | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |  |  |  |
| STARA-ESCD  | Conductivity and DO probe holder for Orion Star<br>A-series armor  |  |  |  |
| STARA-ESPH  | pH electrode holder for Orion Star A-series armor  |  |  |  |
| 1010003   | Optional universal power adapter, 100-240 V, 50/60 Hz  |  |  |  |
| 087010MD  | Orion RDO optical DO probe, 3m cable with stainless steel guard and calibration sleeve                       |  |  |  |
| <b>087020MD</b> Orion RDO optical DO probe, 6m cable with stainle guard and calibration sleeve          |  |  |  |  |
| <b>087030MD</b> Orion RDO optical DO probe, 10m cable with stainless steel guard and calibration sleeve |  |  |  |  |
| <b>087050MD</b> Orion RDO optical DO probe, 15m cable with stainless steel guard and calibration sleeve |  |  |  |  |
| 087100MD  | Orion RDO optical DO probe, 30m cable with stainless steel guard and calibration sleeve                      |  |  |  |
| 083005MD  | Orion polarographic DO probe, 1.5m cable   |  |  |  |
| 080513 Orion polarographic DO probe maintenance kit   |  |  |  |  |
| 080514  | Orion polarographic DO probe electrolyte solution, 60 mL   |  |  |  |
| 080017  | Calibration sleeve for Orion polarographic DO probes   |  |  |  |
| 087003  | Replacement calibration sleeve for Orion RDO optical probes  |  |  |  |
| 087001  | Replacement optical cap for Orion RDO probes   |  |  |  |



рН 📉

ISE \_

ORP \_

Cond

RDO 🗾 DO 🗾





Dissolved Oxygen Portable Meter

Thermo Scientific Orion Star A123 Dissolved Oxygen Portable Meters combine simplicity with accuracy. A large LCD displays dissolved oxygen readings along with temperature. Icons provide quick updates on battery life and completed calibration status. Simple button layout and onscreen messages help with calibration and setup menu choices.

## **Features and Benefits**

- Large, informative screen clearly shows the main information, including DO reading in mg/L or % saturation, temperature in °C or °F, meter mode and battery life
- Automatic temperature, manual salinity input and manual barometric pressure compensation for an accurate result
- Don't miss a reading AUTO-READ™ locks in the stable reading on your screen and ready indicator alerts when readings are stable
- Two calibration options to meet your procedural requirements: water-saturated air method or a Winkler titration
- Accepts Thermo Scientific Orion polarographic dissolved oxygen probes
- Non-volatile memory holds up to 50 data points
- Four AA batteries (included) provide over 2000 hours of operation or purchase the universal power adapter (sold separately) to run meter on AC power
- Perfectly portable, waterproof and protected to take anywhere with a IP67-rated housing
- 3 year meter replacement warranty



# Product Specifications

| Concentration Range 0 to 20 mg/L Resolution 0.01, 0.1 Relative Accuracy ±0.2 % Saturation Range 0 to 200 Resolution 0.1, 1 Relative Accuracy ±2 % Manual Barometric Pressure Correction Manual Salinity Factor Correction Calibration Features Sensor Type Polarographic Range 0 to 50 °C, 32 to 122 °F Resolution 0.1 Temperature Relative Accuracy  ±0.1 Relative Accuracy  ±0.1 Relative Accuracy  ±0.1 Relative Accuracy  ±0.1 Offset Calibration  Yes |                                   |                    | CE   TUV   FCC                        |  |
|--|-----------------------------------|--------------------|---------------------------------------|--|
| Resolution 0.01, 0.1 Relative Accuracy ±0.2 % Saturation Range 0 to 200 Resolution 0.1, 1 Relative Accuracy ±2 %  Manual Barometric Pressure Correction Manual Salinity Factor Correction Calibration Features Sensor Type Polarographic Range 0 to 50 °C, 32 to 122 °F Resolution 0.1 Relative Accuracy Offset Calibration Yes  Water-saturated air, manual (Winkler) Resolution 0.1 Relative Accuracy Offset Calibration Number                          |                                   | ncentration        | ı                                     |  |
| Relative Accuracy #0.2 % Saturation Range 0 to 200 Resolution 0.1, 1  Relative Accuracy #2 %  Manual Barometric Pressure Correction Manual Salinity Factor Correction Calibration Features Sensor Type Polarographic Range 0 to 50 °C, 32 to 122 °F Resolution 0.1  Temperature Relative Accuracy #0.1  Relative Accuracy #0.1  Offset Calibration Vyes  |                                   | Range              | 0 to 20 mg/L                          |  |
| Accuracy   |                                   | Resolution         | 0.01, 0.1                             |  |
| Range  |                                   |                    | ±0.2                                  |  |
| Resolution 0.1, 1  Relative Accuracy ±2 %  Manual Barometric Pressure Correction  Manual Salinity Factor Correction  Calibration Features  Sensor Type Polarographic  Range 0 to 50 °C, 32 to 122 °F  Resolution 0.1  Temperature  Relative Accuracy  Offset Calibration  Offset Calibration  Number   |                                   | Saturation         |                                       |  |
| Dissolved Oxygen  Relative Accuracy  Manual Barometric Pressure Correction  Manual Salinity Factor Correction  Calibration Features  Sensor Type Polarographic  Range 0 to 50 °C, 32 to 122 °F  Resolution 0.1  Temperature  Relative Accuracy  Offset Calibration Offset Calibration Offset Calibration Offset Calibration Ves  |                                   | Range              | 0 to 200                              |  |
| Accuracy Manual Barometric Pressure Correction Manual Salinity Factor Correction Calibration Features Sensor Type Range Resolution Temperature Relative Accuracy Offset Calibration O.1  Relative Accuracy Offset Calibration Offset Calibration O.1  Number  Accuracy Manual 450.0 to 850.0 mm Hg  Vater-saturated air, manual (Winkler)  Date of the calibration O.1  Yes  |                                   | Resolution         | 0.1, 1                                |  |
| Barometric Pressure Correction  Manual Salinity Factor Correction  Calibration Features  Sensor Type Polarographic  Range 0 to 50 °C, 32 to 122 °F  Resolution 0.1  Temperature  Relative Accuracy  Offset Calibration  Number   |                                   |                    | ±2 %                                  |  |
| Salinity Factor Correction  Calibration Features Sensor Type Polarographic  Range 0 to 50 °C, 32 to 122 °F Resolution 0.1  Temperature  Relative Accuracy Offset Calibration  Number   |                                   | rometric<br>essure | 450.0 to 850.0 mm Hg                  |  |
| Features Water-saturated air, manual (Winkler)  Sensor Type Polarographic  Range 0 to 50 °C, 32 to 122 °F  Resolution 0.1  Temperature Relative Accuracy ±0.1  Offset Calibration Ves  |                                   | linity Factor      | 0 to 45                               |  |
| Range  |                                   |                    | Water-saturated air, manual (Winkler) |  |
| Resolution 0.1  Temperature Relative Accuracy ±0.1  Offset Calibration Yes   | :                                 | nsor Type          | Polarographic                         |  |
| Temperature Accuracy ±0.1  Offset Calibration Yes  | 1                                 | inge               | 0 to 50 °C, 32 to 122 °F              |  |
| Accuracy ±0.7  Offset Calibration Yes  | 1                                 | solution           | 0.1                                   |  |
| Calibration Yes  | incialnic                         |                    | ±0.1                                  |  |
| Number   |                                   |                    | Yes                                   |  |
| Datalogging of Datalogs 50   | Datalogging of Datalogs 50        |                    | 50                                    |  |
| <b>Log Function</b> Manual, Automatic with AUTO-READ™  |                                   | g Function         | Manual, Automatic with AUTO-READ™     |  |
| Inputs D0 with Temperature Probe 9-pin mini-DIN  | Inputs Temperature 9-pin mini-DIN |                    | 9-pin mini-DIN                        |  |
| Power Cord Optional – universal, 100-240 VAC   | I                                 | wer Cord           | Optional – universal, 100-240 VAC     |  |
| Power Battery Power Included – 4 AAs   | er I                              | ttery Power        | Included – 4 AAs                      |  |
| Battery Life 2000 hrs  | I                                 | ttery Life         | 2000 hrs                              |  |

## Meter Ordering Information

| Cat. No.  | Description   |
|-----------|---|
| STARA1230 | Orion Star A123 Dissolved Oxygen Portable Meter - Four AA batteries (installed)   |
| STARA1235 | Orion Star A123 Dissolved Oxygen Portable Meter Kit  - 083005MD Orion polarographic DO probe  - 080017 callibration sleeve for DO probe  - 080513 DO probe maintenance kit  - Protective armor with electrode holder  - Hard carrying case  - Four AA batteries (installed) |

## Accessories

| Cat. No.   | Description  |  |  |  |
|------------|--|--|--|--|
| STARA-CS   | Hard carrying case for Orion Star A-series portable meters   |  |  |  |
| STARA-AR   | Armor for Orion Star A-series portable meters, includes electrode holders for pH, conductivity and DO probes |  |  |  |
| STARA-ESCD | Conductivity and DO probe holder for Orion Star<br>A-series armor  |  |  |  |
| STARA-ESPH | pH electrode holder for Orion Star A-series armor  |  |  |  |
| 083005MD   | Orion polarographic DO probe, 1.5 m cable  |  |  |  |
| 080513     | Orion polarographic DO probe maintenance kit   |  |  |  |
| 080017     | Calibration sleeve for Orion polarographic DO probes   |  |  |  |
| 1010003    | Optional universal power adapter, 100-240 V, 50/60 Hz  |  |  |  |



ISE

ORF

Cond

RD(

DO





# **Thermo Scientific Orion Meter Accessories**

For Orion Star A and VERSA STAR Series

With the new Orion meters, no details were too small. Even our accessories have undergone a redesign.

#### Accessories

| Accessories  |  |  |  |  |
|--|--|--|--|--|
| Cat. No.   | Description  |  |  |  |
| STARA-CS   | Hard field case  |  |  |  |
| 1210005  | Soft field case  |  |  |  |
| STARA-AR   | Armor for Orion Star A portable meters (includes pH electrode storage sleeve and holders for pH, conductivity and DO probes) |  |  |  |
| STARA-ESPH   | Replacement pH electrode storage sleeve holder (for use with armor)  |  |  |  |
| STARA-ESCD   | Replacement conductivity and DO storage sleeves holder (for use with armor)  |  |  |  |
| STARA-BEA  | Orion Star A electrode arm and holder head   |  |  |  |
| STARA-HB Weighted, detached base for electrode arm and holder head                           |  |  |  |  |
| STARA-ELHD   | Orion Star A electrode holder head (attaches to existing Orion Star electrode arm)   |  |  |  |
| 1010003  | Universal power adapter (for Orion Star A benchtop and portable meters)  |  |  |  |
| STARA-PWR  | Replacement universal power adapter for<br>Orion VERSA STAR  |  |  |  |
| 1010004 RS232 computer and printer cable set   |  |  |  |  |
| 1010006 Star series inkjet printer, 110V/220V and cable                                      |  |  |  |  |
| 1010008  | 010008 Ink cartridge for Orion Star series printer   |  |  |  |
| 1010009  | 010009 Paper for Orion Star series printer   |  |  |  |
| 1010053 RS232 computer cable   |  |  |  |  |
| Orion Star stirrer probe (automatic stirrer probe wit for benchtop Orion Star series meters) |  |  |  |  |
| 096021   | Replacement paddles for 096019 - 3 pack  |  |  |  |
| 090045   | 990045 BNC shorting cap  |  |  |  |
|  |  |  |  |  |



Armor (cat. no. STARA-AR) adds extra meter protection. The strap can attach to a belt or pipe to keep the meter close to the tester or test site. The built-in stand in the armor works great for additional hands-free operation.



The redesigned electrode holder head (cat. no. STARA-ELHD, sold separately or included with cat. no. STARA-BEA) reduces the amount of sample needed for testing since the probe heights are adjustable.







# **Thermo Scientific Orion Meter Accessories**

For Orion Star, Star Plus, Aplus and PerpHecT Meter Series

To make sure you have everything you need to attain accurate results for your applications, we offer a multitude of accessories for both field and lab applications. From printers to field cases, electrode stands and stirrers, we have it all. Below are the top selling accessories — for complete details, please visit

## www.thermoscientific.com/water.

#### Aplus™ and PerpHecT® Meter Accessories

| Apius    | and respined i weter Accessories   |  |  |
|----------|--|--|--|
| Cat. No. | Description  |  |  |
| 020125   | 110 V power adapter  |  |  |
| 020130   | 220 V power adapter  |  |  |
| 020135   | 240 V power adapter  |  |  |
| 090043   | Swing arm electrode stand  |  |  |
| 090070   | Basic electrode stand and holder   |  |  |
| 910002   | Universal electrode holder for basic electrode stand   |  |  |
| 090045   | BNC shorting cap   |  |  |
| 0900A0   | 900A thermal printer for Aplus meters, 110 V   |  |  |
| 0900A1   | 900A thermal printer for Aplus meters, 220 V   |  |  |
| 0900A7   | 900A thermal printer for Aplus meters, 240 V   |  |  |
| 0900A4   | Paper for 900A printer   |  |  |
| 0900A6   | Battery pack for 900A printer  |  |  |
| PRT300   | PRT300 ink printer and cable for Aplus and PerpHecT meters, 110 V  |  |  |
| PRT301   | PRT300 ink printer and cable for Aplus and PerpHecT meters, 220 V  |  |  |
| PRT302   | Ink ribbon for PRT300 printer  |  |  |
| 0ACBL0   | RS232 computer cable for Aplus and PerpHecT meters   |  |  |
| OPBLSK   | Aplus portable meter starter kit – includes field case, field stand, beaker, 2 x 60 mL bottles each of pH 4.01, 7.00, and 10.01 buffers, and 1 x 60 mL bottle of pH storage solution         |  |  |
| OBCHSK   | Aplus benchtop meter starter kit – includes basic electrode stand and holder, 917005 epoxy ATC probe, and 1 x 475 mL bottle each of pH 4.01, 7.00, and 10.01 buffers and pH storage solution |  |  |
| 0607CA   | 607 multi-channel electrode switchbox and cables, 110 V  |  |  |
| 0607CI   | 607 multi-channel electrode switchbox and cables, 220 V  |  |  |
|          |  |  |  |







Soft Field Case, Cat. No. 1210005

Hard Field Case, Cat. No. 1210004

#### Orion Star™ Series Meter Accessories

| Cat. No.  | Description  |  |  |
|---|--|--|--|
| 1010003   | Universal power adapter  |  |  |
| 1010004   | RS232 computer and printer cable set   |  |  |
| 1010005   | RS232 to USB cable adapter   |  |  |
| 1010006   | Orion Star series inkjet printer, 110 V/220 V and cable  |  |  |
| 1010008   | Ink cartridge for Orion Star series printer  |  |  |
| 1010009   | Paper for Orion Star series printer  |  |  |
| 1010053   | RS232 computer cable   |  |  |
| 1210001   | Portable star series electrode stand   |  |  |
| 1210002   | Meter impact jacket and electrode storage chamber  |  |  |
| 1210003   | Electrode storage chamber (for meter impact jacket)  |  |  |
| 1210004   | Hard field case  |  |  |
| 1210005   | Soft field case  |  |  |
| 096019  | Orion Star stirrer probe (automatic stirrer probe with paddle for benchtop Orion Star series meters)   |  |  |
| 096021 Replacement paddles for 096019 - 3 pack  |  |  |  |
| 090043 Swing arm electrode stand  |  |  |  |
| 090044 Electrode holder head with 15mm dia cond cell slot, for use with 090043                            |  |  |  |
| 090045  | BNC shorting cap   |  |  |
| 1010050 Adapter – MiniDIN ATC probe to 3.5 mm phono tip<br>(310, 320, 330, 350, 370, 555A, 550A, 535A, an |  |  |  |
| 1010051   | Adapter - MiniDIN ATC probe to 8 pin DIN meters (210Aplus, 230Aplus, 250Aplus, 290Aplus, 410Aplus, 420Aplus, 525Aplus, 710Aplus, 720Aplus, and 920Aplus) |  |  |
| 1010800 Adapter – MiniDIN DO probe to 8 pin waterproof DIN Emeters (830A, 835A and 862A)                  |  |  |  |
| 1010801   | Adapter – MiniDIN DO probe to 8 pin DIN DO meters (805Aplus, 810Aplus, and 850Aplus)   |  |  |
| 1010900   | Adapter – MiniDIN conductivity probe to 8 pin waterproof DIN conductivity meters (1230, 555A, 550A, 162A, 135A, and 130A)                                |  |  |
| 1010901   | Adapter – MiniDIN conductivity probe to 8 pin DIN conductivity meters (150Aplus, 145Aplus, 125Aplus, 115Aplus, and 105Aplus)                             |  |  |



# Thermo Scientific Orion Intrinsically Safe\* Meters



## **Features and Benefits**

- For use where non-sparking equipment is desired, these carry US intrinsically safe approval rating by Factory Mutual (FM)\*
- Waterproof, portable models for basic field measurements
- 266S and 136S meters include RS232 output, advanced GLPdoc 100 point datalogger with time and date stamp and includes software program for simplified data collection and remote operation on a computer
- Accepts Waterproof DIN, dual banana plug connections
- Units use 3 AA batteries
- 131S and 136S meters offer linear and non-linear temperature compensation for natural and ultrapure water applications
- 3 year meter warranty
- \* Does not meet European Intrinsically Safe specification.

#### Accessories

| Cat. No. | Description    |  |
|----------|----------------|--|
| 013655   | Electrode clip |  |



## **Product Specifications**

|          |                            | 266S/261S  | 136S/131S                                       |
|----------|----------------------------|--|---|
|          | Range                      | -2.00 to 16.00   | -   |
|          | Resolution                 | 0.1, 0.01  | -   |
|          | Relative<br>Accuracy       | ±0.01  | -   |
| рН       | Auto-Buffer<br>Recognition | Dependent on<br>buffer set (7 pre-<br>programmed sets) | -   |
|          | Calibration<br>Points      | 3  | -   |
|          | Range                      | ±1300.0  | -   |
| mV       | Resolution                 | ±0.1   | -   |
|          | Relative<br>Accuracy       | ±0.1 % of value<br>+0.3 mV                             | -   |
|          | Range                      | -  | 0.00 to 500 mS                                  |
|          | Resolution                 | -  | 0.01 up to 99.99 µS/cm<br>0.1 above 99.99 µS/cm |
| Cond     | Relative<br>Accuracy       | -  | ±0.5 % of measured value + 1 count              |
|          | Calibration<br>Points      | -  | 1   |
|          | Cell Constants             | -  | 0.010 to 199.0 cm <sup>-1</sup>                 |
| Salinity | Range                      | _  | 0 to 45.0<br>(10 to 30 °C) g/kg                 |
|          | Resolution                 | -  | 0.1   |
|          | Range                      | _  | 0 to 1999 mg/L<br>(10 to 40 °C)                 |
|          | Resolution                 | -  | 1   |
|          | Range                      | -20 to 120 °C  | -20 to 120 °C                                   |
| Temp     | Resolution                 | 0.1  | 0.1   |
| Temp     | Relative<br>Accuracy       | ±0.3 °C  | ±0.3 °C   |

## Portable Meter and Electrode

| Cat. No. | Meter | Sensor (s)   |   |
|----------|-------|--|---|
| 0261S0   | 261S  | 9109WP low maintenance pH/ATC Triode $^{\text{\tiny M}}$ |   |
| 0266S0   | 266S  | 9109WP low maintenance pH/ATC Triode                     |   |
| 0131S0   | 131S  | 013010A DuraProbe™ conductivity cell                     | 4 |
| 0136S0   | 136S  | 013010A DuraProbe conductivity cell                      |   |

#### Portable Meter. Electrode and Field Kit

| Fortable Metel, Electrode and Field Kit |       |  |   |  |  |
|---|-------|--|---|--|--|
| Cat. No.                                | Meter | Sensor (s)                                 | Accessories & Solutions   |  |  |
| 0261S3                                  | 261S  | 9109WP low<br>maintenance<br>pH/ATC Triode | Neck strap, field case,<br>stand, electrode holder,<br>beakers, 4, 7 and 10 buffers                                 |  |  |
| 0266\$3                                 | 266S  | 9109WP low<br>maintenance<br>pH/ATC Triode | Neck strap, field case,<br>stand, electrode holder,<br>beakers, 4, 7 and 10<br>buffers, RS232 cable<br>and software |  |  |
| 0131\$3                                 | 131S  | 013010A DuraProbe conductivity cell        | Beakers, standards  |  |  |
| 0136S3                                  | 136S  | 013010A DuraProbe conductivity cell        | Beakers, standards,<br>RS232 cable and software   |  |  |

We have more to present in our newly revamped electrodes & accessories sections.
Keep going!

рН

ISE

ORP

Cond